

THE ONLY HIGH-LEVEL DISINFECTANT PROVEN EFFECTIVE AGAINST HPV IN 30 SECONDS, DESIGNED TO TREAT THE ENTIRE ULTRASOUND STATION

There are over 200 different Human Papillomavirus (HPV) types. Types 16 and 18 are known to be responsible for approximately 70% of cervical cancer cases. A scientific study by Meyers, et al. (2020) demonstrates that Tristel Duo ULT possesses efficacy against infectious HPV types 16 and 18 on a transvaginal ultrasound probe in 30 seconds.

Tristel Duo ULT provides high-level disinfection for endocavity ultrasound probes, skin surface ultrasound probes (including those used during invasive procedures), probe holders, cords, plugs, keyboards, monitors and even mobile medical devices that combine an ultrasound probe with a smart device.



TRISTEL DUO ULT IS EASY TO USE. SIMPLY APPLY, WIPE AND WAIT!



DISPENSE TWO ALIQUOTS OF TRISTEL DUO ONTO A DRY WIPE (DUO WIPE IS RECOMMENDED).



USE THE WIPE TO SPREAD THE FOAM OVER THE SURFACE OF THE ULTRASOUND PROBE AND ASSOCIATED PARTS. ENSURE ALL AREAS ARE COVERED.



LEAVE THE SURFACE TO DRY TO ENSURE A MINIMUM CONTACT TIME OF 30 SECONDS.

Tristel Duo[™]
ULT



EFFECTIVE

Tristel Duo ULT is sporicidal, mycobactericidal, virucidal, fungicidal and bactericidal with a contact time of 30 seconds. The product has been extensively tested according to European Standard suspension tests and a standardised methodology that involves the inoculation of surfaces with the test organism. Tristel's chlorine dioxide is effective against microorganisms of concern such as:

- *Bacillus subtilis*
- *Mycobacterium terrae*
- **Human Papillomavirus (HPV)**
- Hepatitis B Virus (HBV)
- Human Immunodeficiency Virus (HIV)
- Herpes Simplex Virus Type 1
- *Aspergillus brasiliensis*
- *Candida albicans*
- *Escherichia coli*
- *Neisseria gonorrhoeae*
- *Staphylococcus aureus*
- *Gardnerella vaginalis*
- *Streptococcus agalactiae*
- Carbapenem-resistant Enterobacteriaceae (CRE)
- *Klebsiella pneumoniae*
- Vancomycin-resistant Enterococci (VRE)
- *Enterococcus faecium*
- Methicillin-resistant *Staphylococcus aureus* (MRSA)



SIMPLE

Training and certification are an essential part of Tristel's customer service. Our Online Training Portal provides Tristel Duo ULT user training at a time and place convenient to you.



ECONOMICAL

Tristel Duo ULT does not require investment in expensive equipment or service contracts. Two aliquots of foam and one Duo Wipe is all it takes to achieve disinfection to the highest level.



MOBILE

Tristel Duo ULT can be used wherever required, without having to worry about the limitations that automated systems bring. Tristel Duo ULT is also suitable for use in remote locations where electricity and water are not readily available.



COMPATIBLE

Tristel Duo ULT has been tested and confirmed to be compatible by major manufacturers, including:

- Alpinion Medical
- BD (Bard Access)
- BK Medical
- Butterfly Network
- Canon Medical Systems
- Carestream
- Esaote
- Exact Imaging
- FUJIFILM Healthcare
- GE Healthcare
- Healcerion
- MCube
- Mindray
- MobileODT
- Philips
- Quantel Medical
- Samsung Healthcare
- Siemens Healthineers
- Sonoscape
- Supersonic Imagine
- Verathon

APPLICATIONS


Tristel Duo ULT is the only high-level disinfectant that treats all parts of an ultrasound probe the same. It is the only disinfectant that can travel alongside mobile medical devices as an easy-to-use high-level disinfection method.

MONITOR


HOLDER

PLUG

CABLE



TABLET SCREEN



SMART PHONES



HANDLE
SHAFT

INVASIVE
ULTRASOUND PROBES

SKIN SURFACE
ULTRASOUND PROBES



TRISTEL DUO ULT

- 125ml Tristel Base Solution (citric acid) + 125ml Tristel Activator Solution (sodium chlorite)
- When mixed upon pressing the foam pump, Tristel's proprietary chlorine dioxide chemistry is generated.
- 310 aliquots per bottle/2 bottles per box

TRISTEL CLEAN

- 310 aliquots per bottle/6 bottles per box
- 310 aliquots per bottle/2 bottles per box

TRISTEL DUO ULT

DESTROYS HPV ON REAL LIFE DEVICES*

FIND OUT MORE: BIT.LY/HPVARTICLE

*Meyers, C., Milici, J., Robison, R. (2020) "The Ability of Two Chlorine Dioxide Chemistries to Inactivate Human Papillomavirus-contaminated Endocavitary Ultrasound Probes and Nasendoscopes" DOI: 10.1002/jmv.25666.

Duo Wipes are designed specifically for the application of Tristel Duo ULT.

Duo Wipes are made of 100% non-woven, low-linting polypropylene (18gsm).

The use of Duo Wipes allows the effective release of the disinfectant onto the surface.



Scan the QR code to find out more about digital traceability.

Tristel™

UK Manufacturer: Tristel Solutions Limited, Unit 1b, Lynx Business Park, Fortham Road, Snailwell, Cambridgeshire, CB8 7NY, UK
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

India: Tristel India Private Ltd.
T +91 88619 23023 - E mail-india@tristel.com - W www.tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Contact us on
Whatsapp



Join us on
LinkedIn

