


Tristel Fuse™

for Surfaces

High-level disinfecting & sporicidal solution for preventing outbreaks and ward closures




- ☐ Tighten the spigot from the inside before using TANK.
- ☐ Pour 5 litres of ambient water into TANK.

Note: Do not use hot water.

- ☐ Fold the sachet in half and squeeze one side to burst contents through centre seal. Contents will start to turn yellow. Fuse is now ready to be diluted.

- ☐ Tear or cut the sachet. Pour the sachet contents into the water in TANK. Empty the sachet below the water line would help to avoid splashing the concentrate.



- ☐ Screw close the TANK lid, ensuring it is completely closed. The solution can be kept for up to **24 hours**.

- ☐ Record the batch number, date and time of preparation and sign on the label provided.

- ☐ Use the tap on the front of TANK to decant Fuse working solution into smaller volume bottles or containers.

- ☐ Apply Fuse to surfaces with a cloth for damp dusting.
- ☐ Leave surfaces to dry naturally to allow for at least a 5-minute contact time.

Notes:

- Check use by date before use.
- Make up solution in a well-ventilated area.
- Wear appropriate PPE before commencing and wash hands after removing gloves and apron.
- Avoid contact with skin and eyes. Contact with the disinfectant may cause mild irritation. If contact with skin, wash affected areas with plenty of soap and water. In case of contact with eyes, rinse immediately with water and seek medical advice.
- Store in original packaging in a cool, well-ventilated area out of direct sunlight. Keep between 10°C and 35°C.
- Do not use on soft furnishings or fabrics.
- Dispose of empty packaging in accordance with local policy
- TANK should be placed on a non-slip surface and remain stationary when in use.
- Rinse and dry TANK and dosing bottles between uses.
- Do not use if packaging is damaged.
- Do not place objects on top of TANK when in use.
- For professional use only.
- Wipe the labels clean to record new data (for fresh solution preparation).

Tristel

Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 • E mail@tristel.com • W www.tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Hong Kong & Taiwan: 21st Floor, 168 Electric Road, Fortress Hill, Hong Kong
T +852 2895 6968 • F +852 2869 4388 • E customerservicehk@tristel.com

Copyright © Tristel Solutions
Mkt-Wal-1335-1