

reusable hygienic tourniquet



1. application

Loop daisygrip around the patient's limb, click magnetic clasp together.



2. tightening

Place finger between strap and skin to prevent pinching. Pull strap from close to the clasp at a 45° angle to tighten – do not pull from the distal (tail end) of the strap.



3. line insertion/blood collection

Loosen the daisygrip once the needle is in the vein.



4. loosening

Place finger under clasp. Use a 'lifting' motion to flip the clasp in the direction of the strap to loosen daisygrip and release stasis.





5. removing

Loosen, then slide the top clasp sideways with your thumb, whilst holding the bottom arm of the clasp steady (with fingers or palm).



6. reprocessing

Ensure clasp and strap are cleaned and/or disinfected in line with local requirements after every use. Strap should be dry before next use.



Distributed by:

Tristel™

Australian Sponsor:

Tristel Pty Ltd, Unit 44, 328 Reserve Road,
Cheltenham, VIC 3192, Australia
T 1300 680 898 - E mail-au@tristel.com
W www.tristel.com/au-en/

Tristel New Zealand:

23 Birch Avenue, Judea,
Tauranga P.O. Box 3110
T +64 (0) 7 5771560 - E mail-nz@tristel.com
W www.tristel.com/nz-en/

 Always follow the directions for use.

What do you think of daisygrip?

