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TEMPS DE CONTACT DE 30 SECONDES

- NHS. 2018. Ultrasound Scan. [online] Available at: <https://www.nhs.uk/conditions/ultrasound-scan/> [Accessed 3 October 2020].
- Centers for Disease Control and Prevention (2008) Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008 Available at: <https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines-H.pdf> [Accessed 3 October 2020].
- Abramowicz JS, Evans DH, Fowlkes JB, Marsal K, Tierhaar G, on behalf of the WFUMB Safety Committee. Guidelines for cleaning transvaginal ultrasound transducers between patients. *Ultrasound in Med & Biol.* 2017; 43(5):1076-1079.
- EFUMB (European Federation of Societies for Ultrasound in Medicine and Biology). -2017; <http://www.efumb.org/safety/resources/2017-probe_cleaning.pdf>
- Prevention Du Risque Infectieux Associe Aux Actes D'échographie Endocavitaire <https://solidarites-sante.gouv.fr/IMG/pdf/theses_sondes.pdf>
- Najisen, C., Humphreys, H., Koerner, R., Grenier, N., Brady, A., Sidhu, P., Nicolau, C., Mostbeck, G., D'Onofrio, M., Gangi, A. and Claudon, M., 2017. Infection prevention and control in ultrasound - best practice recommendations from the European Society of Radiology Ultrasound Working Group. *Insights into Imaging*, 8(6), pp.523-535.
- Guidelines for reprocessing ultrasound transducers AJUM, 2017, 201:30-40.
- BMUS / SOR (British Medical Ultrasound Society / Society and College of Radiology). -2017; <https://www.bmus.org/policies-statements-guidelines/professional-guidance/guidelines-for-professional-ultrasound-practice/>
- Irish Health Service Executive (HSE) Quality Improvement Division - Decontamination Safety Programme (2017) HSE guidance for decontamination of semi-critical ultrasound probes. Semi-invasive and Non-invasive Ultrasound Probes QPSD-GL-026-1 <http://www.hse.ie/en/g/about/Who/OD/national-safety-programmes/decontamination/Ultrasound-Probe-Decontamination-Guidance-Feb-17.pdf>
- Health Facilities Scotland Decontamination Services (2016) NHS Scotland guidance for decontamination of semi-critical ultrasound probes; semi-invasive and non-invasive ultrasound probes <http://www.hfscotland.nhs.uk/documents/hal/infectioncontrol/guidelines/NHSScotland-Guidance-for-Decontamination-of-Semi-Critical-Ultrasound-Probes.pdf>
- Guidelines for Reprocessing Ultrasound Transducers by the Australasian Society for Ultrasound in Medicine and the Australasian College for Infection Prevention and Control (2017) AJUM 20 (1) <http://online.library.wiley.com/doi/10.1002/ajum.12042.pdf>
- Rutala WA, Weber D. Reprocessing semicritical items. *Am J Infect Control.* 2016;44:e53-e62. doi: 10.1016/j.ajic.2015.12.029.
- Hamm, R., Combs, C. and Davidson, C., 2020. Society for Maternal-Fetal Medicine Special Statement: Reducing the risk of transmitting infection by transvaginal ultrasound examination. *American Journal of Obstetrics and Gynecology*, 223(3), pp.62-66.
- Y'Zai, F., Boumaria, C., Leroy, S., Mekki, Y., Quentin-Noury, C., Kann, M. (2014) Persistence of Microbial Contamination on Transvaginal Ultrasound Probes despite Low-Level Disinfection Procedure. *PLoS ONE*, vol. 9, no. 4 [Online]. DOI: 10.1371/journal.pone.0093368 (Accessed 17 June 2018).
- Ma, S., Yeung, A., Chan, P., Graham, C. (2014) High level disinfection reduces HPV contamination of transvaginal sonography probes in the emergency department [Online]. Available at: <https://emj.bmj.com/content/2016/472.responses/high-level-disinfection-reduces-hpv-contamination-of-transvaginal-sonography-probes-in-the-emergency-department> (Accessed 26 July 2018).
- Casallegro, J.S., Canval, K., Eibach, D., Valdeyron, M.L., Lamblin, G., Jacquemoud, H., Mellier, G., Lina, B., Gaucherand, P., Mathevet, P., Mekki, Y. (2012) High Risk HPV Contamination of Endocavity Vaginal Ultrasound Probes: An Underestimated Route of Nosocomial Infection? *PLoS ONE*, vol. 7, no. 10 [Online]. DOI: doi:10.1371/journal.pone.0048197 (Accessed 23 June 2018).
- Strauss, S., Sastry, P., Sonnex, C., Edwards, S., Gray, J. (2002) Contamination of environmental surfaces by genital human papillomaviruses¹.
- Roden, R., Lowy, D., Schiller, J. (1997) 'Papillomavirus Is Resistant to Desiccation'. *The Journal of Infectious Diseases*, vol. 176, no. 5, pp. 1076-1079 [Online]. DOI: <https://doi.org/10.1086/516515>
- Basseal, J., Westerway, S. and Hjelt, J., 2020. Analysis of the integrity of ultrasound probe covers used for transvaginal examinations. *Infection, Disease & Health*, 25(2), pp.77-81.
- Westerway, S. C., Basseal, J. M., Brockway, A., Hjelt, J. A., Carter, D. A. (2016) 'Potential Infection Control Risks Associated with Ultrasound Equipment – A Bacterial Perspective'. *Ultrasound in Medicine & Biology* [Online]. DOI: 10.1016/j.ultrasmedbio.2016.09.004 (Accessed 25 October 2018).
- Kiegs, M., Sim, B., Thom, G., Tunbridge, M., Barnett, A., Fraser, J. (2019) 'Efforts to Attenuate the Spread of Infection (EASI): a prospective, observational multicentre survey of ultrasound equipment in Australian emergency departments and intensive care units'. *Critical care and resuscitation*, vol. 17, no. 1, pp. 43-46 [Online]. Available at: <https://www.ncbi.nlm.nih.gov/pubmed/25702761> (Accessed 04 April 2019).
- Allos, B.M., Blaser, M.J., Platts-Mills, J. and Kosek, M., 2014. *Campylobacter Species - Infectious Disease and Antimicrobial Agents*. [online] Antimicrobe.org. Available at: <http://www.antimicrobe.org/new/b91.asp> [Accessed 5 August 2020].
- Meyers, J., Ryndock, E., Conway, M., Meyers, C. and Robison, R., 2014. Susceptibility of high-risk human papillomavirus type 16 to clinical disinfectants. *Journal of Antimicrobial Chemotherapy*, 69(6), pp.1546-1550.
- Meyers, C., Milici, J., Robison, R. (2020) 'The Ability of Two Chlorine Dioxide Chemistries to Inactivate Human Papillomavirus-contaminated Endocavitary Ultrasound Probes and Nasendoscopes'. *Published IN THE JOURNAL OF MEDICAL VIROLOGY* [Online]. Available at: <bit.ly/HPVARTICLE>
- World Health Organization. 2018. Human Papillomavirus (HPV). [online] Available at: <https://www.who.int/immunization/diseases/hpv/en/> [Accessed 19 October 2020].