**Tristel Duo** TM/MC

---

**OPH**

**High-level disinfectant foam for semi-critical  
medical devices used in ophthalmology.**

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## INTENDED USE

Tristel Duo OPH is a high-level disinfectant foam intended to disinfect cleaned, reusable, non-lumened, semi-critical medical devices used in ophthalmology. Intended for use by qualified healthcare personnel.

## INDICATIONS FOR USE

Tristel Duo OPH is a chemical disinfectant intended for the high-level disinfection of cleaned, re-usable, non-lumened, semi-critical medical devices used in ophthalmology such as diagnostic lenses, tonometer prisms, pachymeters, A-scan and B-scan biometry probes.

Tristel Duo OPH is mycobactericidal, virucidal, fungicidal and bactericidal with a disinfection contact time of 2 minutes at 20°C and when used according to the Directions for Use.

Tristel Duo OPH must be applied on the surface of a medical device using a Duo Wipe and rinsed from the device using sterile water or a Tristel Rinse Wipe.

## ACTIVE INGREDIENT


Tristel Duo OPH is made of two solutions. Tristel Duo Activator solution and Tristel Duo Base solution are mixed when dispensed to generate the active ingredient **chlorine dioxide** ( $\text{ClO}_2$ ). The in-use concentration of chlorine dioxide is 0.032% (320 ppm). The maximum concentration of chlorine dioxide after a 2-minute contact time is 0.14% (1,400 ppm). Tristel Duo Activator and Base solutions are contained in a dual compartment bottle with a dispenser pump and are ready to use.

### Tristel Duo Activator

Key Ingredient:	
Sodium chlorite.....	0.5%
Other ingredients.....	99.5%
Total.....	100%

### Tristel Duo Base

Key Ingredient:	
Citric acid.....	5.26%
Other ingredients.....	94.74%
Total.....	100%

## FIRST AID TREATMENT

Tristel Duo Base contains citric acid.

Do not get in eyes. Wear eye protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice.

KEEP OUT OF REACH OF CHILDREN.

Tristel Duo OPH has not been tested for dermal penetration/absorption or for systemic toxicity resulting from dermal absorption.

## CONTRAINDICATIONS

- Do not use if the Duo bottle or Duo Wipes pouch is damaged, for longer than the prescribed contact time, or on non-compatible devices.
- Tristel Duo OPH efficacy may be compromised if it is used with wipes other than Duo Wipes, and if user instructions are not followed in full.
- Not intended for single use devices.
- Tristel Duo OPH should not be used to sterilize medical devices.

## PRECAUTIONS

Follow OSHA Bloodborne Pathogens Universal Precautions when handling and cleaning soiled devices.

- Wear gloves when reprocessing devices.
- Do not use Tristel Duo OPH with the following materials: aluminium, anodized aluminium, brass, copper, chromium, rubber (NBFQ) or polyurethane elastomer 35A. Testing has shown that Tristel Duo OPH is not compatible with these materials.
- Do not use Tristel Duo OPH beyond 6 months after opening.

## DIRECTIONS FOR USE

### CLEANING


Blood, other bodily fluids, gels, and lubricants must be thoroughly cleaned from the surfaces of medical devices before high-level disinfecting with Tristel Duo OPH. Refer to the reusable device manufacturer's labelling and instructions for procedures on disassembly, handling, cleaning, and decontamination.

### CLEANING AGENT COMPATIBILITY

Tristel Duo OPH is compatible with enzymatic detergents, low foaming, and easily rinsed from the devices. Tristel Clean is recommended.

### OPENING AND SETTING UP THE DUO BOTTLE

- Wear gloves.
- Depress the pump 4 times to prime and prepare for use. Dispensed foam can be flushed with water to drain or contained on a wipe and disposed to waste in accordance with local regulations. Tristel Duo OPH is now ready for use.
- Use Tristel Duo OPH within 6 months after opening.
- Record the date the bottle is first used, i.e. opened, and the 'Use by date' in the designated area on the front label of the bottle.
- Check that the 'Use by date' does not extend past the expiration date. The expiration date is printed on the Tristel Duo OPH label. If the 6 months from opening extends beyond the expiration date, record the expiration date in the 'Use by date' field on the label.

## HIGH-LEVEL DISINFECTION

### STEP 1

- Check the 'Use by date' on the Duo bottle. Do not use if out of date.
- Remember, if the Duo bottle is being used for the first time, the pump must be primed. Refer to the instructions above for opening and setting up the Duo bottle.
- Always use Duo on a cleaned device.
- Wear new gloves when handling Tristel Duo OPH and medical devices.

### STEP 2

Dispense 3 aliquots of Tristel Duo onto a Duo Wipe or directly onto the device. Use the wipe to spread the foam over the surface of the device and ensure all areas are covered. Continuous wiping of the surface is not necessary and should be avoided to prevent shedding of the wipe.

**Note:** Discard the used wipe in accordance with local regulations. Do not macerate. Do not re-use.

### STEP 3

After wiping, place the device on a clean surface to avoid recontamination. The contact time for high-level disinfection is **2 minutes**.

It is recommended that a timer is used to ensure the contact time of 2 minutes is observed.

### STEP 4

Rinse the device by wiping with a Tristel Rinse Wipe or immersing it completely in sterile water. One rinse is required to remove Tristel Duo OPH residues unless more rinse cycles are specified by the device manufacturer. Always follow device manufacturer's labelling for rinsing instructions.

Thoroughly dry all surfaces of the device using a clean soft, lint-free wipe or cloth.

## **STORAGE CONDITIONS AND EXPIRATION DATE**

Store Tristel Duo OPH in their original containers at temperature 15-30°C in a ventilated area out of direct sunlight.

Once opened, use Tristel Duo OPH within 6 months, ensuring that this does not extend past the expiration date. The expiration date is indicated on the front label.

## **EMERGENCY AND ADDITIONAL INFORMATION**

For further safety information please refer to the Safety Data Sheet. Emergency information can be obtained from Tristel in the UK at + 44 (0)1638 721500.

## **DISPOSAL INFORMATION**

Do not reuse or refill Duo bottle.

Dispose of the empty container in accordance with municipal/provincial/territorial requirements. Offer for recycling, if available.

## PRODUCT OPTIONS

DESCRIPTION	CASE CONTENTS	ORDER CODE
Tristel Duo OPH	2 Duo Bottles 1 User Guide	TSL023925
Tristel Duo OPH	6 Duo Bottles 1 User Guide	TSL022725

## ADDITIONAL PRODUCTS

DESCRIPTION	CASE CONTENTS	ORDER CODE
Duo Wipes - each pouch contains 160 dry wipes	2 Pouches	TSL032401
Duo Wipes - each pouch contains 160 dry wipes	6 Pouches	TSL032301
Tristel Clean - each bottle contains 240ml enzymatic cleaning solution, delivered in ready-to-use foam	6 Bottles 1 User Guide	TSL023301
Tristel Rinse Wipes - containing 50 x pre-wetted, sterile rinse wipes	50 Sachets 1 User Guide	TSL030301
Duo DOK - storage solution for Tristel Duo OPH	5 Units 1 User Guide	TSL023401

**Tristel Duo OPH is classified as a Class II Medical Device under the Canadian Medical Device Regulations SOR/98-282.**

**Duo Wipes, Tristel Clean and Tristel Rinse Wipes are classified as Class I Medical Devices under the Canadian Medical Device Regulations SOR/98-282.**

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Tristel Solutions Ltd, Lynx Business Park,
Tristel Solutions Limited, Lynx Business Park,