
3

Tristel™

Rinse Wipe

USER INSTRUCTIONS

- Disinfect hands and wear gloves when handling disinfectants and medical devices.
- Wear appropriate Personal Protective Equipment (PPE).
- Store between 10°C and 35°C in the original packaging. Do not freeze.
- Do not use if the Wipe is damaged, discoloured, or dry.
- Keep the empty Wipe sachet(s) to complete the Tristel Quality Audit Trail Record Book at the end of the decontamination procedure.
- Continuous wiping of the surface is not necessary and should be avoided to prevent shedding of the Wipe.
- For professional use only.
- These Instructions for Use (IFU) should be used in conjunction with the product label, Safety Data Sheet (SDS) and device manufacturer's instructions.
- Do not use if the Wipe sachet is damaged.
- Do not use past the expiry date given on the sachet. Discard expired Wipes according to local regulations.
- For batch number and expiry date see sachet.

Intended Use

The third and final step in the decontamination procedure is rinsing of the medical device. The Rinse Wipe is impregnated with de-ionised water and a low-level of antioxidant which removes chemical residues from a surface. Each Rinse Wipe sachet is packed and then sterilised by gamma-irradiation.

STEP 1

Take one Rinse Wipe sachet.

STEP 2

Remove the Wipe from its sachet and lay it out in the palm of your hand.

STEP 3

Wipe the surface of the device that has been decontaminated to remove excess foam.

STEP 4

Discard the used Wipe and gloves in accordance with local regulations. Do not reuse the Wipe. Keep the empty Wipe sachet for traceability.

Note: Upon completion of the decontamination cycle the device should be left to air dry. Store the device in accordance with hospital protocols to prevent damage or recontamination.

Please complete the Tristel Quality Audit Trail Record Book to keep record of the decontamination procedure.

If you believe that this device has contributed to a serious incident (i.e. a death, the temporary or permanent serious deterioration of a person's state of health, or a serious public health threat), contact the device manufacturer and your national competent authority immediately.

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IMPORTERS

**EUROPEAN UNION (EU) / EUROPEAN ECONOMIC AREA (EEA)
/ NORTHERN IRELAND / SWITZERLAND / TURKEY**

**AUSTRIA, BULGARIA, CROATIA, CZECH REPUBLIC,
ESTONIA, GERMANY, GREECE, HUNGARY, LATVIA,
LITHUANIA, POLAND, SLOVAKIA, SLOVENIA**

Tristel GmbH

Karl-Marx-Allee 90A, 10243 Berlin, Germany
T +49 (0)30 54844226 E deutschland@tristel.com
W www.tristel.de

BELGIUM, FRANCE, LUXEMBOURG, NETHERLANDS

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NORTHERN IRELAND

BRENNAN & CO NI

Unit 6 & 7, Carrowreagh Business Park, Dundonald, Belfast, Northern Ireland, BT16 1QQ
T +44 2890 410 136 E info@brennanconi.com
W www.brennanconi.com

NORWAY

VINGMED AS

Solbråveien 13, N-1383 Asker, Norway
T +47 67 58 06 80 E info@vingmed-as.no
W www.vingmed-as.no

SPAIN

VESISMIN, S.L.

C/ Aribau 230-240, 6ª planta, L-O, 08006 Barcelona, Spain
T +34 934 09 53 01 E vesismin@vesismin.com
W www.vesismin.com

SWITZERLAND

TRISTEL AG

Sandgrube 29, 9050 Appenzell, Switzerland
T +41 (0)71 5670658 E schweiz@tristel.com
W www.tristel.com

TURKEY

YEFE DIS TICARET LTD STI.

Cennetoğlu Mah. Yeşillik Cad. 228 Karabağlar, İzmir 35400, Turkey
T +90 532 275 1133 E efeayvali@pelisLtdcom
W www.yefe.com.tr

For any further information, please contact us at mail@tristel.com.

Tristel™

WE HAVE CHEMISTRY.

Manufactured by: Tristel Solutions Limited
Lynx Business Park, Fordham Road, Snailwell, Cambridgeshire CB8 7NY
T +44 (0)1638 721500 - F +44 (0)1638 721911 - E mail@tristel.com - W www.tristel.com

Australian sponsor: Tristel Pty Ltd. 40/328 Reserve Road, Cheltenham, VIC 3192, Australia
T 1300 680 898 - E mail-au@tristel.com

Tristel Rinse Wipes is entered into the Australian Register of Therapeutic Goods as a Class I medical device, ARTG number: 173046.

For Tristel patent information please visit www.our-patents.info/tristel

Contact Tristel, your local distributor or visit www.tristel.com for supporting documents such as safety data sheets and reports.

Where necessary, external translations have been sourced from ISO 17100:2015 approved supplier: www.tolingo.com.

EC REP Tristel NV, Smallandlaan 14B,
2660 Antwerp, Belgium

CH REP Tristel AG, Sandgrube 29,
9050 Appenzell, Switzerland

RIN/098/5
14-JUN-2022

STERILE R