Tristel Fuse™

for Medical Surfaces

High-level disinfecting & sporicidal solution

High-level disinfecting and sporicidal solution for medical surfaces

APPLICATIONS AND USES

- Medical surfaces within operating theatres
- Commodes
- Incubators
- IV poles
- Patient trolleys
- Hospital mattresses and bed frames
- Other hard non-porous medical surfaces including floors, walls, doors

Note: Tristel Fuse for Medical Surfaces is not suitable for use on invasive medical devices.

EACH SACHET CONTAINS

- 50ml Base solution (5% solution of citric acid in demineralised water)
- 50ml Activator solution (2.1% sodium chlorite solution in demineralised water)
- This makes 100ml chlorine dioxide in aqueous solution for dilution into 5 litres of water.

THE TRISTEL CHEMISTRY

Tristel Fuse for Medical Surfaces utilises Tristel's patented chlorine dioxide chemistry, a well-documented and highly effective biocide. The chemical symbol for chlorine dioxide is ClO₂.

MODE OF ACTION

Tristel Fuse for Medical Surfaces incorporates two separate compartments that contain the Tristel base and activator solutions. When mixed by bursting the sachet, chlorine dioxide is generated. Chlorine dioxide is a powerful oxidising agent – an electron receiver. This means that the chlorine dioxide molecule is in constant search for an additional electron. When a bacterial cell comes into contact with chlorine dioxide, it donates an electron from its cell wall. This creates a breach in the cell wall through which cell contents pass in an attempt to bring the concentrations on either side of the cell membrane to equilibrium. The cell dies through lysis.

BIOCIDAL PERFORMANCE

Tristel Fuse for Medical Surfaces is sporicidal, mycobactericidal, virucidal, fungicidal and bactericidal with a contact time of only 5 minutes. Tristel Fuse for Medical Surfaces is effective against microorganisms of concern such as:

Spores	<i>Clostridium difficile</i> , <i>Bacillus cereus</i> , <i>Bacillus subtilis</i>
Mycobacteria	<i>Mycobacterium avium</i> , <i>Mycobacterium terrae</i>
Viruses	Adenovirus Type 5, Poliovirus Type 1, Murine Norovirus
Fungi	<i>Aspergillus brasiliensis</i> , <i>Candida albicans</i>
Bacteria	<i>Enterococcus hirae</i> , <i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i> , Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA), <i>Escherichia coli</i> , Multi-drug resistant <i>Acinetobacter baumannii</i> (MDRAB), Vancomycin-resistant Enterococci (VRE) <i>Enterococcus faecium</i> , Carbapenem-resistant Enterobacteriaceae (CRE) <i>Klebsiella pneumoniae</i>

EN Activator

Warning

- Contact with acids liberates very toxic gas.
- Causes skin irritation.
- Causes serious eye irritation.

Base Warning


- Causes skin irritation.
- Causes serious eye irritation.
- Harmful to aquatic life with long lasting effects.

Activator and Base Precautionary statements:

- Wear protective gloves/protective clothing/eye protection.
- IF ON SKIN:** Wash with plenty of water/.
- IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- If skin irritation occurs:** Get medical advice.
- If eye irritation persists:** Get medical advice.


HOW TO USE TRISTEL FUSE FOR MEDICAL SURFACES

- Make up solution in a well ventilated area.
- Ensure appropriate PPE according to hospital protocol is worn prior to preparing Fuse.
- Store in original packaging in a cool, well ventilated area out of direct sunlight.
- Once mixed with water, solution can be kept for up to 24 hours when stored in an air tight container.
- Dispose of empty packaging in accordance with local policy and national regulations.
- Avoid contact with skin and eyes. Contact with the disinfectant may cause mild irritation. Wash affected areas with plenty of soap and water.
- For professional use only.
- Wear an apron and nitrile gloves before commencing and wash hands after removing apron and gloves.

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

Note: Do not use hot water.


Take one sachet.

Fold sachet in half and squeeze one side to burst contents through centre seal. Contents will start to turn yellow.

Tear or cut sachet.

Take care when opening sachet, do not splash concentrate. Keep away from face and eyes.

Product is ready to use.

Apply to surface with mop for floors or a cloth for damp dusting other surfaces and leave to dry naturally for a 5 minute contact time.

Manufactured by: Tristel Solutions Ltd

For Tristel patent information please visit www.our-patents.info/tristel

Contact Tristel, your local distributor or visit www.tristel.com for supporting documents such as safety data sheets, microbiological test data and reports.

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Wyprodukowano przez: Tristel Solutions Ltd.
Tristel Solutions Ltd, Lynx Business Park, Snailwell, Cambridgeshire, CB8 7NY, UK