
Tristel Duo™

ORL

Tristel Duo ORL is a chlorine dioxide based high-level disinfectant foam. It is CE marked as a Class IIb Medical Device (MDD 93/42 EEC).

INTENDED USE

Tristel Duo ORL is designed specifically for the high-level disinfection of non-lumened invasive medical devices used in otorhinolaryngology (ORL).

THE TRISTEL CHEMISTRY

Tristel Duo ORL is based on Tristel's proprietary chlorine dioxide (ClO₂) chemistry, a well-documented and highly effective biocide. Tristel Duo ORL incorporates two separate compartments that contain Tristel Base Solution (citric acid) and Tristel Activator Solution (sodium chlorite). When mixed upon pressing the foam pump, Tristel's proprietary chlorine dioxide chemistry is generated. Tristel Duo ORL contains a non-ionic surfactant for added cleaning performance.

BIOCIDAL PERFORMANCE

Tristel Duo ORL is sporicidal, mycobactericidal, virucidal, fungicidal and bactericidal with a contact time of 30 seconds. The product has been tested according to the relevant European Standards required for disinfectants used in medical areas. Tristel Duo ORL is effective against microorganisms of concern such as:

- *Bacillus subtilis*
- *Clostridium difficile*
- *Mycobacterium terrae*
- Human Papillomavirus (HPV)
- Human immunodeficiency virus (HIV)
- Hepatitis B virus (HBV)
- *Vaccinia virus*
- *Candida albicans*
- *Aspergillus brasiliensis*
- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*
- *Enterococcus hirae*
- Carbapenem-resistant Enterobacteriaceae (CRE) *Klebsiella pneumoniae*
- Vancomycin-resistant Enterococci (VRE) *Enterococcus faecium*

USER INSTRUCTIONS

- Do not use if the bottle is damaged.
- Disinfect hands and wear appropriate Personal Protective Equipment (PPE) when handling disinfectants and medical devices.
- If Tristel Duo ORL is used for the first time, depress the pump two to four times to prime it.
- For use on pre-cleaned devices only.
- For professional use.
- These Instructions for Use (IFU) should be used in conjunction with the product label, Safety Data Sheet (SDS) and device manufacturer's instructions.

Step 1

Dispense two aliquots of Tristel Duo onto a dry wipe (Duo Wipe is recommended). Use the wipe to spread the foam over the surface of the medical device and ensure all areas are covered. Continuous wiping of the surface is not necessary and should be avoided to prevent shedding of the wipe.

Note: Discard the used wipe in accordance with local regulations. Do not macerate. Do not re-use.


Step 2

Leave the surface to dry to ensure a minimum contact time of 30 seconds.

Step 3

Rinse the instrument with water of appropriate quality.

Note: Store the medical device in accordance with hospital protocols to prevent damage or recontamination.

Warning

Base

Causes serious eye irritation.

Precautionary statements

Base

Wear eye protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If eye irritation persists:
Get medical advice.

Manufactured by: Tristel Solutions Ltd

For Tristel patent information please visit www.our-patents.info/tristel

Contact Tristel, your local distributor or visit www.tristel.com for supporting documents such as safety data sheets, microbiological test data and reports.

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Geproduceerd door: Tristel Solutions Ltd
Tristel B.V., Binderij 7 R, 1185 Amstelveen, Nederland
Tristel Solutions Ltd, Lynx Business Park, Fordham Road, Snailwell, CB8 7NY, UK