# Tristel Fuse™

## for Surfaces

High-level disinfecting & sporicidal solution

# High-level disinfecting and sporicidal solution for surfaces

## APPLICATIONS AND USES

- Tristel Fuse for Surfaces is intended for disinfection of hard non-porous surfaces.

## EACH SACHET CONTAINS

- 50ml Base solution (5% solution of citric acid in demineralised water)
- 50ml Activator solution (2.1% sodium chlorite solution in demineralised water)
- This makes 100ml chlorine dioxide in aqueous solution for dilution with water.

## THE TRISTEL CHEMISTRY

Tristel Fuse for Surfaces utilises Tristel's patented chlorine dioxide chemistry, a well-documented and highly effective biocide. The chemical symbol for chlorine dioxide is ClO<sub>2</sub>.

## MODE OF ACTION

Tristel Fuse for Surfaces incorporates two separate compartments that contain the Tristel Base and Activator solutions. When mixed by bursting the sachet, chlorine dioxide is generated. Chlorine dioxide is a powerful oxidising agent – an electron receiver. This means that the chlorine dioxide molecule is in constant search for an additional electron. When a bacterial cell comes into contact with chlorine dioxide, it donates an electron from its cell wall. This creates a breach in the cell wall through which cell contents pass in an attempt to bring the concentrations on either side of the cell membrane to equilibrium. The cell dies through lysis.

## BIOCIDAL PERFORMANCE


Tristel Fuse for Surfaces is sporicidal, mycobactericidal, virucidal, fungicidal and bactericidal with a contact time of five minutes.

Tristel Fuse for Surfaces is effective against microorganisms of concern such as:

Spores	<i>Clostridium difficile</i> <i>Bacillus cereus</i> <i>Bacillus subtilis</i>
Mycobacteria	<i>Mycobacterium avium</i> <i>Mycobacterium terrae</i>
Viruses	Adenovirus Type 5 Poliovirus Type 1 Murine Norovirus
Fungi	<i>Aspergillus brasiliensis</i> <i>Candida albicans</i> <i>Candida auris</i>
Bacteria	<i>Enterococcus hirae</i> <i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Streptococcus pyogenes</i> Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) <i>Escherichia coli</i> Multi-drug resistant <i>Acinetobacter baumannii</i> (MDRAB) Vancomycin-resistant Enterococci (VRE) <i>Enterococcus faecium</i> Carbapenem-resistant Enterobacteriaceae (CRE) <i>Klebsiella pneumoniae</i>


## HOW TO USE TRISTEL FUSE FOR SURFACES

- Make up solution in a well ventilated area.
- Ensure appropriate PPE according to hospital protocol is worn prior to preparing Fuse.
- Wear gloves and apron before commencing and wash hands after removing gloves and apron.
- Store in original packaging in a cool, well ventilated area out of direct sunlight.
- Once mixed, solution can be kept for up to 24 hours when stored in an air tight container.
- Dispose of empty packaging in accordance with local policy and national regulations.
- Avoid contact with skin and eyes. Contact with the disinfectant may cause mild irritation. Wash affected areas with plenty of soap and water.
- For professional use only.

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

*Note: Do not use hot water.*

**2** Take one sachet.


**3** Fold sachet in half and squeeze one side to burst contents through centre seal. Contents will start to turn yellow.

**4** Product is ready to use. Tear or cut sachet.

Take care when opening sachet, do not splash concentrate. Keep away from face and eyes.

**5** Pour contents into the five litres of ambient water.

**6** Apply to surface with mop for floors or a cloth for damp dusting other surfaces and leave to dry naturally.

**Manufactured by:** Tristel Solutions Ltd

For Tristel patent information please visit [www.our-patents.info/tristel](http://www.our-patents.info/tristel)

Contact Tristel, your local distributor or visit [www.tristel.com](http://www.tristel.com) for supporting documents such as safety data sheets, microbiological test data and reports.

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Tillverkad av: Tristel Solutions Ltd.
Levererad av: Vingmed AB, Box 576, SE-175 26 JÄRFÄLLA, Sweden
Tristel Solutions Ltd, Lynx Business Park, Snailwell, Cambridgeshire, CB8 7NY, UK