




High-level disinfecting & sporicidal solution for preventing outbreaks and ward closures.

Tristel Fuse for Surfaces combines the powerful biocidal efficacy of chlorine dioxide with cleaning action. Designed specifically for use on large surface areas such as floors and walls.

Effective against a wide range of microorganisms, including SARS-CoV-2, the virus which causes COVID-19, and multi-drug resistant pathogens, in **five minutes**.

EFFECTIVE AGAINST SARS-CoV-2



Pour five litres of ambient water into a container.
Note: Do not use hot water.



Take one sachet and fold it in half. Squeeze one side to burst contents through the centre seal. Contents will start to turn yellow.



Tear or cut the sachet corner and pour contents into five litres of water. Tristel fuse for surfaces is now ready to use.



Apply tristel fuse for surfaces with a mop for floors and walls and with a cloth for damp dusting. Observe a contact time of **five minutes**.

Tristel FuseTM
for Surfaces

Classified as biocide under the EU Biocides Regulation. Use biocides safely. Always read the label and product information before use.

WHY CHOOSE TRISTEL FUSE FOR SURFACES?

Tristel Fuse for Surfaces combines the powerful biocidal efficacy of chlorine dioxide with cleaning action. It is designed specifically for use on most floors, walls and high-grade stainless steel surfaces.

Tristel Fuse for Surfaces is delivered in a unique dual-compartment burstable sachet containing Tristel Base and Activator solutions. When mixed upon bursting the sachet, Tristel's proprietary chlorine dioxide chemistry is generated.

Each sachet produces five litres of working solution at **one concentration**, with one contact time to destroy a wide range of microorganisms. Tristel Fuse for Surfaces achieves high-level disinfection, including sporicidal efficacy, in **five minutes**.

Training and certification are an essential part of Tristel's customer service. Our Online Training Portal provides Tristel Fuse for Surfaces user training at a time and place convenient to you.

Effective against microorganisms of concern such as:

- *Bacillus subtilis*
- Mycobacteria causing Tuberculosis (TB)
- SARS-CoV-2
- Poliovirus
- Adenovirus
- Norovirus
- *Candida albicans*
- *Candida auris*
- *Escherichia coli*
- *Pseudomonas aeruginosa*
- *Enterococcus hirae*
- *Staphylococcus aureus* (including Methicillin-resistant *Staphylococcus aureus* (MRSA))
- Carbapenemase-resistant Enterobacteriaceae (CRE) *Klebsiella pneumoniae*
- Vancomycin-resistant Enterococci (VRE) *Enterococcus faecium*
- Multi-drug resistant *Acinetobacter baumannii* (MDRAB)

CHLORINE DIOXIDE

Tristel Fuse for Surfaces utilises Tristel's proprietary chlorine dioxide chemistry (ClO₂), a well-documented and highly effective biocide. ClO₂ is a strong oxidant whose germicidal characteristics are well known.


Literature indicates that it causes the chemical disruption of cell walls, damage to inner cell membranes, the denaturing of proteins and the impairment of genetic material.

Tristel Fuse for Surfaces is superior to surface cleaners and disinfectants based on bleach, chlorine tablets, quaternary ammonium compounds, phenols and aldehydes.

APPLICATIONS

Tristel Fuse for Surfaces is designed specifically for the cleaning and high-level disinfection of large hard non-porous surfaces in critical areas, such as:

- Operating Theatres
- Oncology Units
- Burn Units
- Intensive Care Units
- General Wards
- Laboratories
- Care Environment
- Accident & Emergency
- Emergency Vehicles



TRISTEL FUSE FOR SURFACES

- 40 sachets/box
- Order code: **TSL010201**
- Classified as biocide under the EU Biocides Regulation. Use biocides safely. Always read the label and product information before use.

Tristel™

UK Manufacturer: Tristel Solutions Limited,
Unit 1b, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

Hong Kong: Tristel Asia Limited, 21st Floor, 168 Electric Road,
Fortress Hill, Hong Kong - T +852 2895 6968 - F +852 2869 4388
E customerservicehk@tristel.com

Copyright
© Tristel Solutions
Mkt-Bro-T546-3
14/MAY/2024