



Semi-automated disinfection made simple

Stella is a semi-automated washer disinfectant designed specifically for the disinfection of invasive and non-invasive, heat sensitive, non-lumened and single-lumened medical devices used in Urology, Gynaecology, IVF, Cardiology, ENT, Anesthesiology and Respiratory.

Stella combines the simplicity of manual soaking with the sophistication of a fully automated washer disinfectant.

APPLICATIONS

High-level disinfection of semi-critical medical devices, including:

- Bronchoscopes
- Cystoscopes
- Hysteroscopes
- Intubation endoscopes
- Laryngoscopes
- Manometry catheters
- Nasendoscopes
- Transesophageal Echocardiography (TOE/TEE) probes
- Transrectal probes
- Transvaginal probes

COMPATIBLE

Tristel Fuse for Stella has been tested and confirmed to be compatible with medical devices from manufacturers, including*:

- BK Medical
- Canon Medical Systems
- Carestream
- FUJIFILM Healthcare
- Karl Storz
- Laborie
- Philips
- Samsung Healthcare

*Material compatibility refers to certain medical devices/models.

PUBLICATIONS

- ☰ Stella reduces capital spend and reprocessing times. Xuan Wu Hospital, China. January 2016.
- ☰ Stella enables Shanghai Pulmonary Hospital to high-level disinfect bronchoscopes nine times faster. Shanghai Pulmonary Hospital, China. January 2016.
- ☰ A randomised single-blind comparison of the effectiveness of the high-level disinfectants Tristel Fuse (chlorine dioxide) and Cidex OPA (ortho-phthalaldehyde) for use with flexible cystoscopes. Peter J Gilling, Michael Addidle, Rana Reuther, Michelle Lockhart, Christopher Frampton and Mark R Fraundorfer. 2013.
- ☰ Stella for simple cystoscope reprocessing. Time for a change at the Reading Urology Partnership. April 2012.
- ☰ Stella for simple GI manometry disinfection. Rising to the challenge at Castle Hill. May 2012.



INCREASES PATIENT THROUGHPUT

Stella's high-level disinfection cycle time is five minutes. When cleaning is included, the cycle time is ten minutes.



MOBILE & LIGHTWEIGHT

Stella measures 70cm x 48cm x 18cm and weighs 5.5kg when empty. The base and lid can be stacked and easily transported due to their portable size and weight.



GIVES YOU PEACE OF MIND

Prior to the start of each cycle, Stella confirms that Tristel solutions for Stella are used. It also checks that there is no channel blockage.

At the end of each cycle, Stella provides the user with a unique traceability code.



ECONOMICAL

Stella requires minimal capital investment when compared to a fully automated washer disinfectant. The short cycle time improves the turnaround of medical devices enabling you to operate with fewer of them.



SAFE FOR MEDICAL DEVICES

When Stella irrigates the lumen, it monitors the back pressure of the Tristel solution passing through it and adjusts the pressure to protect it. The solution also pulses through the lumen, rather than moving at an even flow, dislodging soil and bio burden in the process.



IMPROVES HEALTH & SAFETY

Tristel Fuse for Stella utilises Tristel's proprietary chlorine dioxide (ClO₂) chemistry, a well-documented and highly effective biocide. The Health and Safety profile of chlorine dioxide is superior to that of sodium hypochlorite, peracetic acid and glutaraldehyde.




Refer to the Stella Help Guides, Quick Guides and Instructions for Use for full instructions.

KEY

Stella System A

Stella System B

Stella System C



1 Turn on Stella IQ.

2 Turn on Stella Pulse.



3 Add the cleaned medical device – connect lumened medical devices (in the case of Stella B/C) to Stella Pulse.

Note: When using Stella C, the medical device should be pre-cleaned.



4 If using Stella for medical device cleaning prior to high-level disinfection, prepare five litres of Tristel Clean for Stella and add to the inner compartment of the base.



5 For the high-level disinfection cycle, add one sachet of Tristel Fuse for Stella to five litres of water and add to the inner compartment of the base.



6 Close the lid and Stella will do the rest. The high-level disinfection cycle time is five minutes. When cleaning is included, the overall cycle time is ten minutes. Stella drains automatically.



7 At the end of the cycle, Stella will issue a validation code to confirm it has completed successfully. This should be recorded in the Stella Quality Audit Trail Record Book or downloaded to the Stella Suite.

PRODUCT INFORMATION

STELLA SYSTEM A

For the high-level disinfection of non-lumened medical devices.



STELLA SYSTEM B

For the high-level disinfection of non-lumened and single-lumened medical devices.

STELLA SYSTEM C

For the cleaning and high-level disinfection of non-lumened and single-lumened medical devices.

More product information, such as safety data sheets, reports, publications on materials and studies, is available on request or online at www.tristel.com

TRISTEL FUSE FOR STELLA

40 sachets per box



TRISTEL CLEAN FOR STELLA

1 litre cleaning solution

Tristel™

UK Manufacturer: Tristel Solutions Limited, Unit 1b, Lynx Business Park, Fordham Road, Snailwell, Cambridgeshire, CB8 7NY, UK
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

Australia: Tristel Pty Ltd, 44/328 Reserve Road, Cheltenham, VIC 3192
T 1300 880 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

Hong Kong: Tristel Asia Ltd, 21st Floor, 168 Electric Road, Fortress Hill, Hong Kong
T +852 2895 6968 - F +852 2869 4388 - E customerservicehk@tristel.com

Singapore: Tristel Pte Limited, 2 Venture Drive, Vision Exchange #14-04, Singapore 608526
T +65 6518 9918 - E mail-sg@tristel.com

New Zealand: Tristel New Zealand Limited, 23 Birch Avenue, Judea, Tauranga
T +64 (0)7 5771560 - F +64 (0)7 5771567 - E mail-nz@tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>
Copyright © Tristel Solutions - Mkt-Fac-791 - 29/MAY/2024 - V8