

# reusable hygienic tourniquet



## 1. application

Loop daisygrip around the patient's limb, click magnetic clasp together.



## 2. tightening

Place finger between strap and skin to prevent pinching. Pull strap from close to the clasp at a 45° angle to tighten – do not pull from the distal (tail end) of the strap.



## 3. line insertion/blood collection

Loosen the daisygrip once the needle is in the vein.



## 4. loosening

Place finger under clasp. Use a 'lifting' motion to flip the clasp in the direction of the strap to loosen daisygrip and release stasis.



## 5. removing

Loosen, then slide the top clasp sideways with your thumb, whilst holding the bottom arm of the clasp steady (with fingers or palm).



## 6. reprocessing

Ensure clasp and strap are cleaned and/or disinfected in line with local requirements after every use. Strap should be dry before next use.

Distributed by:

**Tristel**™

Tristel Solutions Ltd

Unit 1b, Lynx Business Park, Fordham Road, Snailwell, Cambridgeshire, CB8 7NY, UK  
T +44 (0) 1638 721500 – E mail@tristel.com – W www.tristel.com

Copyright © Tristel Solutions Limited – Mkt-Wal-2794 – 31/JUL/2024 – V1



Always follow the directions for use.

What do you think of daisygrip?

