Tristel Trio[™]

Wipes System

Quality Audit Trail Record Book

THE TRISTEL TRIO WIPES SYSTEM

The Tristel Trio Wipes System is a three-part system for high-level disinfection of non-lumened invasive and non-invasive medical devices. It comprises three wipes and an Activator Foam that in sequence perform the steps of the decontamination procedure.

1 Tristel Pre-Clean Wipe

Impregnated with enzymatic detergent to clean instruments as the first part of the decontamination process.

2a Tristel Sporicidal Wipe 2b Tristel Activator Foam

A high-level disinfectant wipe which generates chlorine dioxide when activated with Tristel Activator Foam.

3 Tristel Rinse Wipe

Impregnated with de-ionised water and a low-level of antioxidant to rinse instruments as the final stage of decontamination process. The wipe is sterile packed.

Quality Audit Trail Record Book

Finally, the Quality Audit Trail Record Book records the decontamination process through to use on a patient or to storage of the device.

It tracks the following steps:

- Identification of the medical device to be decontaminated and the date and time it is decontaminated.
- Cleaning of the medical device (with the Pre-Clean Wipe or by other means).
- Correct activation and use of the Sporicidal Wipe.
- Rinsing of the medical device (with the Rinse Wipe or by other means).
- Identification of the next destination for the device – patient or storage.
- Identification of the person responsible for the decontamination process.


USING THE QUALITY AUDIT TRAIL RECORD BOOK

The Record Book will be maintained in the work area where the device(s) will be decontaminated.

All decontamination procedures will be recorded in the Record Book. A new section of the Record Book will be used for each decontamination procedure – i.e. one device per section.

The person decontaminating the device will perform all stages of the decontamination procedure and will then complete the Record Book.

Device being decontaminated

Completion of this section of the Record Book will include:

- Confirm the device being decontaminated.
- A record of its serial number **SN** or other means of identification.
- A record of the date **📅** and time **🕒** of the decontamination procedure.

Using the Pre-Clean Wipe **1**

If the Pre-Clean Wipe has been used, confirm **✓** in the box and log the batch number **LOT** and use by date **📅** as indicated on the sachet used. If the Pre-Clean Wipe has not been used, confirm **✗** in box and state other method employed.

Using the Sporicidal Wipe **2a 2b**

When using the Sporicidal Wipe, write the Activator Lot No. **LOT** onto the **traceability label**, then peel off and affix the 'Record Book Label' part of the **traceability label** to the box identified in the book (2a Sporicidal Wipe).

To activate the Sporicidal Wipe, follow the detailed instructions for use that are printed on the Sporicidal Wipe sachet and also in the user guide. Confirm in the **2b Activator Foam** section that the following steps have taken place:

Remove the Wipe from its sachet and lay it out in the palm of your hand. Remove the lid from the Activator Foam bottle. Apply the required aliquots of Activator Foam onto the Sporicidal Wipe. Confirm **✓** in box.

- If you are using Trio Wipes System **30**, apply two aliquots of Activator Foam onto the Sporicidal Wipe.
- If you are using Trio Wipes System **3**, apply four aliquots of Activator Foam onto the Sporicidal Wipe.

Fold the Wipe in on itself and scrunch together **15 seconds** to activate. Ensure that the Wipe is evenly covered with foam. Confirm **✓** in box.

Wipe the surface of the device in one movement to cover it with foam, ensuring all areas come into contact with the Wipe. Pay special attention to edges, ridges and indentations. Observe a **30-second** contact time. Confirm **✓** in box.

Using the Rinse Wipe **3**

If the Rinse Wipe has been used, confirm **✓** in the box and log the batch number **LOT** and use by date **📅** as indicated on the sachet used. If the Rinse Wipe has not been used, confirm **✗** in box and state the other method employed.

Destination of device **📍**

If the device is to be used on a patient, the 'Patient's Notes Label' part of the **traceability label** should be peeled off and attached to the Patient's notes.

If the device is to be stored and decontaminated again before next use, peel off and affix the 'Patient's Note Label' part of the **traceability label** to the box identified in the book **📍**.

Person Responsible **👤**

The person responsible for performing and completing the decontamination procedure must print their full name **👤** and sign **✍️** the Record Book.

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The Trio Wipes System is classified under the Australian Register of Therapeutic Goods as a Class IIb medical device, AUST R number: 182843.

For Tristel patent information please visit www.our-patents.info/tristel
Contact Tristel, your local distributor or visit www.tristel.com for supporting documents such as safety data sheets, microbiological test data and reports.

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Tristel Solutions Limited, Lynx Business Park, Fordham Road, Snailwell, Cambridgeshire, UK, CB8 7NY
Tristel Pty Ltd. 40/328 Reserve Road, Cheltenham, VIC 3192, Australia