


Enzystel™

Enzystel is a triple-enzymatic detergent for the manual cleaning of endoscopes and surgical instruments including soaking in ultrasonic baths. It removes and dissolves proteinaceous and organic matter such as blood, saliva and mucous from instruments quickly and effectively.



HOW TO USE FOR MANUAL CLEANING




1 Use one 20ml bottle dose in 1 litre of water for instrument soaking (1 part Enzystel in 50 parts water).




2 Soak the instrument for 10 minutes and rinse thoroughly.
Note: Solution should be changed weekly or when it becomes cloudy/soiled.

HOW TO USE FOR AUTOMATED SYSTEMS



1 Use one 10ml bottle dose in 1 litre of water for a flushing solution (1 part Enzystel in 100 parts water).



2 Use for 5 minutes or as per manufacturer's recommended machine times.
Note: Solution should be changed weekly or when it becomes cloudy/soiled.

PRODUCT OPTIONS



Enzystel 1L

MD315

Tristel™
WE HAVE CHEMISTRY.

Manufactured in the United Kingdom by:
Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>