
2a+b

Tristel™

Sporicidal Wipe
& Activator Foam

The second step in the decontamination procedure is high-level disinfection of the medical device. The Sporicidal Wipe is CE Marked as a Class IIb Medical Device (MDD 93/42 EEC).

APPLICATIONS

Applications of the Tristel Trio Wipes System include (but are not limited to) the decontamination of nasendoscopes, transoesophageal echocardiogram (TOE/TEE) probes, transvaginal and transrectal ultrasound probes and other invasive ultrasound probes, non-invasive ultrasound probes (including those used during invasive procedures), laryngoscopes, intubation endoscopes, manometry catheters and ophthalmic devices.

THE TRISTEL CHEMISTRY

The Sporicidal Wipe incorporates Tristel's proprietary chlorine dioxide (ClO₂) chemistry, a well-documented and highly effective biocide. The Wipe is impregnated with Tristel Base Solution (citric acid) and the Tristel Activator Foam is a dilute solution of sodium chlorite. When mixed upon applying Activator Foam onto the Sporicidal Wipe and scrunching them together, Tristel's proprietary chlorine dioxide chemistry is generated.

BIOCIDAL PERFORMANCE OF THE TRISTEL SPORICIDAL WIPE

The Sporicidal Wipe destroys organisms of concern such as bacterial spores, mycobacteria, viruses, fungi and bacteria with a contact time of 30 seconds. It has been extensively tested according to European Standard suspension tests and a standardised methodology that involves the inoculation of surfaces with the test organism. Tristel's chlorine dioxide is effective against microorganisms of concern such as:

- *Clostridium sporogenes*
- Human Papilloma Virus (HPV)
- *Mycobacterium tuberculosis*
- Poliovirus Type 1
- Human Herpesvirus Type 1
- Adenovirus Type 5
- *Candida albicans*
- *Staphylococcus aureus*
- Vancomycin-resistant Enterococci (VRE) *Enterococcus faecium*
- *Klebsiella pneumoniae*
- *Escherichia coli*

USER INSTRUCTIONS

- Disinfect hands and wear gloves when handling disinfectants and medical devices.
- Wear appropriate Personal Protective Equipment (PPE).
- Do not use if the Wipe sachet is damaged or the Wipe is discoloured, damaged or dry.
- Do not use the Activator Foam if the bottle is damaged.
- Keep the empty Wipe sachets to complete the Tristel Quality Audit Trail Record Book at the end of the decontamination procedure.
- Continuous wiping of the surface is not necessary and should be avoided to prevent shedding of the Wipe.
- For professional use only.
- These Instructions for Use (IFU) should be used in conjunction with the product label, Safety Data Sheet (SDS) and device manufacturer's instructions.
- Do not use past the expiry date. For the expiry date and LOT number please see the surface of the carton.

- 1 Disinfect your hands and put on new gloves.
- 2 Take one Sporicidal Wipe sachet.
- 3 Remove the Wipe from its sachet and lay it out in the palm of your hand.

- 4 Remove the lid from the Activator Foam bottle. Apply two aliquots of Activator Foam onto the Sporicidal Wipe.

Note: Activate the Tristel Sporicidal Wipe as soon as you have removed it from the sachet and use it immediately.

Note: If the Activator Foam bottle is being used for the first time, depress the pump two to four times to prime it. The first output from the foam bottle can be left on the Wipe, to be followed by complete aliquots. The Activator Foam bottle is then primed for subsequent use.

- 5 Fold the Wipe in on itself and scrunch together **15 seconds** to activate. Ensure that the Wipe is evenly covered with foam. Presence of a chlorine-like odour confirms that the Wipe is ready to use.

Use the activated Wipe immediately.

- 6 Wipe the surface of the medical device in one movement to cover it with foam, ensuring all areas come into contact with the Wipe. Pay special attention to edges, ridges, indentations, and areas where different materials connect.

- 7 Observe a **30-second** contact time.
- 8 Discard the used Wipe in accordance with local regulations. Do not reuse. Keep the empty Wipe sachet for traceability.

**DO NOT
SHAKE THE
ACTIVATOR
FOAM
BOTTLE.**

Please complete the Tristel Quality Audit Trail Record Book to keep record of the decontamination procedure.

Manufactured by: Tristel Solutions Limited

Lynx Business Park, Fordham Road, Snailwell, Cambridgeshire CB8 7NY

T +44 (0)1638 721500 - F +44 (0)1638 721911 - E mail@tristel.com - W www.tristel.com

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Suisse: Tristel AG, Sandgrube 29, 9050 Appenzell - T +41 (0)71 5670658

E schweiz@tristel.com

A series of horizontal lines for writing, spaced evenly down the page.

A series of 25 horizontal grey lines are spaced evenly down the page, providing a template for writing.

A series of horizontal lines for writing, spaced evenly down the page.

Tristel™

WE HAVE CHEMISTRY.

Tristel Solutions Limited, Lynx Business Park, Fordham Road, Snailwell, Cambridgeshire, CB8 7NY
T +44 (0)1638 721500 - F +44 (0)1638 721911 - E mail@tristel.com - W www.tristel.com

Australian sponsor: Tristel Pty Ltd. 40/328 Reserve Road, Cheltenham, VIC 3192, Australia
T 1300 680 898 - E mail-au@tristel.com

The Tristel Sporidical Wipes & Activator Foam are classified under the Australian Register of Therapeutic Goods as a Class IIb medical device, AUST R number: 182843.

For Tristel patent information please visit www.our-patents.info/tristel

Contact Tristel, your local distributor or visit www.tristel.com for supporting documents such as safety data sheets, microbiological test data and reports.

EC REP Tristel NV, Smallandlaan 14B,
2660 Antwerp, Belgium

CE
2797

TRL/348/5
January 2020

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Fabriqué par: Tristel Solutions Ltd