

Tristel Duo™

OPH

Are your disinfection practices compliant with AS/NZS 5369:2023?

- ✓ Approved for ophthalmic devices
- ✓ Two-minute high-level disinfection
- ✓ Effective against Adenovirus
- ✓ Accessible at point-of-care



Compliant disinfection for tonometers, pachymeters, and other devices that touch the eye.

Why Duo OPH?

Australia and New Zealand Standard require semi-critical devices – such as those that contact the eye – to undergo high-level disinfection with a medical device disinfectant found on the Australian Register of Therapeutic Goods (ARTG). Commonly seen products – such as alcohol wipes, or soaks in dilute bleach – are not recognised as medical device disinfectants, and have gaps in testing, evidence, and non-standardised instructions for use. With the implementation of AS/NZS 5369 in December 2023, small operators are classed as Health Service Organisations, and are required to comply.

The technology behind the function

- Tristel Duo OPH is powered by a proprietary chlorine dioxide (ClO_2) formulation. As an oxidising agent, ClO_2 destroys microorganisms via electron exchange, causing a molecular imbalance and leading to its inactivation and/or cell death.
- Tristel ClO_2 is generated at the point of use; each press of the pump delivers a consistent, measured dose of separated base and activator that combine upon delivery.



How to use

To use, apply 2x pumps of foam to a dry wipe, and apply to the device. Allow to remain in contact for 2-minutes prior to rinsing.



Tristel™

Australia:
Tristel Pty Ltd, Unit 44, 328 Reserve Road,
Cheltenham, VIC 3192, Australia T 1300 680 898
E mail-au@tristel.com W www.tristel.com/au-en/

New Zealand:
23 Birch Avenue, Judea, Tauranga 3110
T +64 (0) 7 5771560 E mail-nz@tristel.com
W www.tristel.com/nz-en/

