


Tristel Duo™

NCU

Intermediate Level Disinfectant



Tristel Duo NCU is a chlorine dioxide based intermediate level disinfectant foam. It is approved as an Instrument Grade Intermediate Level Disinfectant under the Australian Government Department of Health Therapeutic Goods Administration and in accordance with the European Medical Devices Directive 93/42/EEC and the 2007/47/EC amendments thereto.

INTENDED USE

Tristel Duo NCU is designed specifically for the intermediate level disinfection of non-invasive ultrasound probes, cables, plugs, probe holders, monitors and control panels.

Note: *not suitable for use on invasive medical devices.*

THE TRISTEL CHEMISTRY

Tristel Duo NCU is based on Tristel's proprietary chlorine dioxide (ClO_2) chemistry, a well-documented and highly effective biocide. Tristel Duo NCU incorporates two separate compartments that contain Tristel Base Solution (citric acid) and Tristel Activator Solution (sodium chlorite). When mixed upon pressing the foam pump, Tristel's proprietary chlorine dioxide chemistry is generated. Tristel Duo NCU contains a non-ionic surfactant for added cleaning performance.

BIOCIDAL PERFORMANCE

Tristel Duo NCU is efficacious against mycobacteria, viruses, fungi and vegetative bacteria. The product has been tested according to the relevant European Standards required for disinfectants used in medical areas. Tristel Duo NCU is effective against microorganisms of concern such as:

- *Mycobacterium terrae*
- *Mycobacterium avium*
- Human Papillomavirus (HPV)*
- Herpes Simplex Virus Type 1
- Poliovirus Type 1
- Adenovirus Type 5
- Vaccinia virus
- Norovirus
- *Aspergillus brasiliensis*
- *Candida albicans*
- Carbapenem-resistant Enterobacteriaceae (CRE)
Klebsiella pneumoniae
- Vancomycin-resistant Enterococci (VRE)
Enterococcus faecium
- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*
- *Enterococcus hirae*

*Based on Polyomavirus SV40 surrogate testing and live infectious HPV strains 16 and 18 as tested by Meyers et al., (2020).

USER INSTRUCTIONS

- Do not use if the bottle is damaged.
- Disinfect hands and wear appropriate Personal Protective Equipment (PPE) when handling disinfectants and medical devices.
- If Tristel Duo NCU is used for the first time, depress the pump two to four times to prime it.
- Not to be used on soft metals such as brass or copper.
- For professional use.

Step 1

Depress the pump to dispense Tristel Duo onto a dry wipe (Duo Wipe is recommended) or directly onto the surface.

Step 2

Use the wipe to spread the foam over the surface and ensure all areas are wet.

Step 3

Allow the surface to remain wet for 30 seconds. Leave to air dry. Discard the used wipe in accordance with local regulations. Do not macerate. Do not re-use.



Warning

Base

Causes serious eye irritation.

Precautionary statements

Base

Wear eye protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If eye irritation persists:
Get medical advice.



Manufactured by: Tristel Solutions Ltd

For Tristel patent information please visit
www.our-patents.info/tristel

Contact Tristel, your local distributor or visit
www.tristel.com for supporting documents
such as safety data sheets, microbiological
test data and reports.

Tristel™

WE HAVE CHEMISTRY.



Tristel Solutions Ltd, Lynx Business Park, Fordham Road, Snailwell, CB8 7NY, UK
+44 (0)1638 721500 - mail@tristel.com - www.tristel.com



Tristel NV, Smalandlaan 14 B,
2660 Antwerp, Belgium

Australian sponsor: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192, AU
1300 680 898 - mail-au@tristel.com - www.tristel.com

New Zealand: Tristel New Zealand Ltd, 23 Birch Avenue, Tauranga 3110, NZ
+64 (0)7 5771560 - info@tristel.co.nz - www.tristel.com

