

Tristel Jet™ FOAM

SPORICIDAL FOAM

Jet is a powerful sporicidal foam for the cleaning and disinfection of environmental surfaces and non-porous materials in laboratories.


Jet is effective against viruses (including coronaviruses and influenza viruses), bacterial spores, protozoa, mycobacteria, fungi and vegetative bacteria in a 2 minute contact time.

Jet destroys DNA and RNA in a matter of seconds.

Jet contains a non-ionic surfactant for added cleaning performance.

Jet is the best safeguard wherever the risks of infection are highest. Jet is compatible with hard, non-porous materials such as laminate, glass, plastic and stainless steel.

Ideal for use on non-porous high-touch surfaces such as:



**EFFECTIVE
AGAINST
CORONAVIRUSES**

Jet achieves bactericidal, virucidal, fungicidal and mycobactericidal efficacy in 30 seconds, and sporicidal efficacy in 2 minutes, including against:

BACTERIA	CONTACT TIME	VIRUSES	CONTACT TIME
<i>Staphylococcus aureus</i>	30 seconds	Poliovirus Type 1	30 seconds
<i>Enterococcus hirae</i>	30 seconds	Adenovirus Type 5	30 seconds
<i>Pseudomonas aeruginosa</i>	30 seconds	Murine Norovirus	30 seconds
Carbapenem Resistant Enterobacteriaceae (CRE) <i>Klebsiella pneumoniae</i>	30 seconds	Influenza Virus (H1N1)	30 seconds
Vancomycin Resistant Enterococci (VRE) <i>Enterococcus faecium</i>	30 seconds	Coronavirus (SARS-CoV-2)*	30 seconds
		Vaccinia Virus	30 seconds
		Human Papillomavirus (using Polyoma virus SV40 Surrogate)	30 seconds
		Human Papillomavirus**	30 seconds
		Herpes Simplex Virus Type 1	30 seconds
MYCOBACTERIA	CONTACT TIME	SPORES	CONTACT TIME
<i>Mycobacterium terrae</i>	30 seconds	<i>Clostridium sporogenes</i>	2 minutes
<i>Mycobacterium avium</i>	30 seconds		
FUNGI	CONTACT TIME		
<i>Candida albicans</i>	30 seconds		
<i>Aspergillus brasiliensis</i>	30 seconds		

Easy to use with any available wipe, to reduce your plastic use!



3.5 Kilos


=



0.8 Kilos

Saving space, reducing waste

Jet can be used with any sustainably sourced wipe, avoiding the use and waste of plastic wipes.



Powered by chlorine dioxide foam

Jet is built upon the same core chlorine dioxide technology as the Tristel medical device disinfectant products, but specially adapted for surface disinfection where cleaning and biofilm destruction are so important.

Chlorine dioxide is widely regarded as one of the most effective disinfectants for biofilm removal and prevention.

Jet delivers chlorine dioxide as a foam and – unlike sprays – generates no aerosols.

Complete compliance with EN 14885:2018 and the latest efficacy regulatory requirements

SPORES	EN 17126	YEASTS/FUNGI	EN 13624
	EN 13704		EN 14562
MYCOBACTERIA	EN 14348	BACTERIA	EN 16615
	EN 14563		EN 13727
VIRUSES	EN 14476		EN 14561
	EN 16615		EN 16615

AVAILABLE TO ORDER NOW!



Each pack contains: 2x Jet Heads and 2x 800ml Jet Packs

One 800ml Jet pack delivers 570 doses of foam

Quote: TSL051814

Jet is classified as biocide under the EU Biocides Regulation. Use biocides safely. Always read the label and product information before use.