

Anistel

FALCON™



FOR THE RAPID DISINFECTION OF ANIMAL HEALTHCARE ENVIRONMENTS

When used as part of a daily or weekly deep clean regime, FALCON is an excellent preventative method for outbreaks.

- Ready-to-use: no dilution required
- Proprietary chlorine dioxide chemistry
- Sporicidal in one minute
- For small hard surface areas
- No rinsing required
- Effective against coronaviruses*



Classified and labelled in accordance with the Biocidal Products Directive 528/2012

WHY CHOOSE FALCON?

FLEXIBLE

Designed for the disinfection of instruments and equipment used within the veterinary environment, including isolation areas, reception areas and operating theatres.

Note: FALCON is not suitable for use on fabric or soft furnishings

EFFECTIVE

FALCON is sporicidal, mycobactericidal, virucidal, fungicidal and bactericidal in a contact time of one minute. It is effective against microorganisms of concern in animal healthcare, including:

- *Bacillus subtilis*
- *Bacillus pumilus*
- *Mycobacterium avium*
- *Mycobacterium terrae*
- Bovine enterovirus Type 1
- Canine parvovirus
- *Feline calicivirus*
- Adenovirus Type 5
- Poliovirus Type 1
- *Candida albicans*
- *Aspergillus brasiliensis*
- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*
- *Enterococcus hirae*
- *Dermabacter hominis*
- *Staphylococcus capitis*
- *Propionobacter acnes*
- *Micrococcus luteus*
- *Staphylococcus epidermidis*
- *Escherichia coli*
- *Klebsiella pneumoniae*
- *Proteus hauseri* (formerly *Proteus vulgaris*)
- *Bordetella bronchiseptica* (kennel cough)
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Multidrug resistant *Acinetobacter baumannii* (MDRB)
- Vancomycin-resistant Enterococci (VRE) *Enterococcus faecium*

EASY TO USE

STEP ONE

Assemble the FALCON bottle. Take a FALCON trigger head and following the red and blue indicators on the bottle label, insert the tubes into the bottle necks. Fasten the levers on each side of the trigger head to secure. FALCON is ready for use.

STEP TWO

Apply directly on to the surface, or on to a dry wipe.

Three aliquots will cover an area of 0.2m².

STEP THREE


Wipe to disinfect the surface.

STEP FOUR

Ensure all areas remain wet for the one-minute contact time. Leave to air dry.

STEP FIVE

Once the FALCON bottle is empty, unfasten the levers and remove the FALCON trigger head. This trigger head can now be used with the second FALCON bottle.



POWERED BY CHLORINE DIOXIDE

FALCON is built upon the same core chlorine dioxide (ClO₂) technology as the Tristel medical device decontamination products, trusted by the NHS for nearly three decades.

FALCON utilises a reusable trigger head which combines sodium chlorite with an organic acid solution, generating ClO₂ at point of use and at one concentration.

ClO₂ achieves its potent biocidal effect through oxidation. ClO₂ oxidises lipids and proteins present in cell membranes, leading to a loss in membrane integrity and ultimately cell death. ClO₂ is widely regarded as one of the most effective disinfectants for biofilm removal and prevention.



TR117

- 1 x trigger head
- 2 x bottles containing 400ml Tristel Base Solution (citric acid) + 400ml Tristel Activator Solution (sodium chlorite).
- 570 doses per bottle.

BROUGHT TO YOU BY TRISTEL IN THE UNITED KINGDOM

Tristel is a manufacturer of products for infection prevention and contamination control, headquartered in Snailwell, Cambs (UK). Please visit www.tristel.com for more information.