



A mobile colposcope needs a disinfectant that is equally mobile.

Based on Tristel's proprietary chlorine dioxide chemistry, Tristel Duo EVE is designed specifically for the high-level disinfection of non-invasive medical devices such as colposcopes and non-porous surfaces of medical equipment used in Women's Health. It does not require electricity or water to do its job. Tristel Duo EVE is effective against microorganisms of concern with a contact time of only 30 seconds.

TRISTEL DUO EVE IS EASY TO USE. SIMPLY APPLY, WIPE AND WAIT!



**DISPENSE TRISTEL DUO EVE
ONTO A DRY WIPE.**

Duo Wipes are recommended.



**SPREAD THE FOAM OVER THE SURFACE
AND ENSURE ALL AREAS ARE WET.**



**ALLOW THE SURFACE TO REMAIN
WET FOR 30 SECONDS.
LEAVE TO AIR DRY.**

Tristel DuoTM
EVE



EFFECTIVE

Based on Tristel's proprietary chlorine dioxide chemistry, Tristel Duo EVE is effective in short contact times against microorganisms of concern in Women's Health, such as:

- Human Papillomavirus (HPV)
- Hepatitis B virus (HBV)
- Hepatitis C virus (HCV)
- Human immunodeficiency virus (HIV)
- *Escherichia coli*
- *Neisseria gonorrhoeae*
- *Staphylococcus aureus*
- *Candida albicans*
- Herpes simplex virus
- *Gardnerella vaginalis*
- *Streptococcus agalactiae*
- Carbapenem-resistant Enterobacteriaceae (CRE)
- *Klebsiella pneumoniae*
- Vancomycin-resistant Enterococci (VRE)
- *Enterococcus faecium*
- Methicillin-resistant *Staphylococcus aureus* (MRSA)



ECONOMICAL

Tristel Duo EVE does not require investment in expensive equipment or service contracts.



SIMPLE

Tristel Duo EVE is applied onto surfaces with a dry wipe.



COMPATIBLE

Tristel Duo EVE has been tested and confirmed to be compatible by major manufacturers, including:

- Alpinion Medical Systems
- Bard
- BK Ultrasound
- Canon (Toshiba)
- Carestream
- Fujifilm Sonosite
- GE Healthcare
- Hitachi
- Healcerion
- Mindray
- Mobile ODT
- Philips
- Samsung
- Supersonic Imaging

Human Papillomavirus (HPV) is a small, non-enveloped virus that demonstrates resistance to a number of disinfectants¹. Clinical research studies have proven that high and low-risk strains of HPV deoxyribonucleic acid (DNA) and HPV-free virions are found on the surfaces of fomites and medical devices within the clinical environment^{2,3,4,5,6,7}. Devices particularly at risk to contamination are colposcopes, which are used for the examination of the cervix and surrounding tissue. The handle of the device is a high-touch area that has shown to be positive for HPV DNA in 70% of samples^{2,4}.

HPV is highly resistant to extreme desiccation on environmental surfaces and retains infectivity for at least

seven days⁸. If devices are inadequately disinfected and HPV-free virions are present, there is potential of HPV transmission between patients, clinicians and devices. Tristel Duo EVE is proven to destroy HPV⁹.

The presence of DNA on the surfaces of medical devices can interfere with sample collection. This is especially relevant in the case of Sexual Assault Nurse Examiner (SANE) practices. Contamination of genetic material from extraneous sources (e.g., from the examiner or device) is one of the greatest concerns in DNA studies¹⁰. Tristel Duo EVE destroys all traces of DNA and ribonucleic acid (RNA).

References

¹ DOI: 10.1093/jac/dku006

³ DOI: [https://doi.org/10.1016/0020-7292\(90\)90371-Q](https://doi.org/10.1016/0020-7292(90)90371-Q)

⁴ DOI: 10.1136/sextrans-2014-051977

² DOI: 10.1136/sti.78.2.135

⁵ DOI: 10.1371/journal.pone.0093368

⁶ DOI: 10.1136/emered-2012-201407

⁷ DOI: doi:10.1371/journal.pone.0048137

⁸ DOI: <https://doi.org/10.1086/516515>

⁹ RVV DKI polyomavirus SV40 testing

¹⁰ DOI: <http://dx.doi.org/10.1155/2015/365674>

Duo Wipes are designed specifically for the application of Tristel Duo EVE.

Duo Wipes are made of 100% non-woven, low-linting polypropylene (18gsm) and pre-cut to prevent shedding when dispensed from the tub.

The use of Duo Wipes allows the effective release of the disinfectant onto the surface.



DUO WIPES

- 200 wipes per tub / 6 tubs per box

TRISTEL DUO EVE

- 125ml Tristel Base Solution (citric acid) + 125ml Tristel Activator Solution (sodium chlorite)

When mixed upon pressing the foam pump, Tristel's proprietary chlorine dioxide chemistry is generated.

- 310 aliquots per bottle / 6 bottles per box

Tristel[™]
WE HAVE CHEMISTRY.

Manufactured in the United Kingdom by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

Hong Kong & Taiwan: 21st Floor, 168 Electric Road, Fortress Hill, Hong Kong
T +852 2895 6968 - F +852 2869 4388 - E customerservice@tristel.com.hk

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Copyright © Tristel Solutions
Mkt-Bro-564-4
07/SEP/2022