



## A CHLORINE DIOXIDE (ClO<sub>2</sub>) SPORICIDAL FOAM FOR THE DISINFECTION OF NON-POROUS MEDICAL SURFACES AND NON-INVASIVE MEDICAL DEVICES.

Not intended to be used with ultrasound probes or invasive medical devices.

### HOW TO USE



1 To assemble JET, take a trigger, unscrew the caps on the bottles, and following the red A and blue B indicators on the bottle label, insert the tubes into the respective bottle necks.



2 Fasten the levers anti-clockwise on each side of your trigger to secure. Ensure the nozzle is aligned to the top line on the trigger head. JET is now ready for use.



3 Apply at least 2 doses onto a wipe.



4 Wipe onto the surface.



5 Ensure the surface remains wet for the 1 minute contact time. Leave to air dry.



6 When you can't recycle, reuse! Read the 'BIG IMPACT' section below to find out more.

### PRE-USE CHECKS

Some important reminders before starting:

- Always read the label, safety data sheets (SDS) and product information before use.
- Wear the appropriate personal protective equipment (PPE).
- Do not use if product is damaged.
- JET is for professional use only.
- Do not apply to brass, copper, aluminium, mild steel, or silver surfaces.
- Carefully clean, rinse and dry the surface before use.



Refer to the User Guide for full instructions.

### BIG IMPACT

Your JET trigger should be reused with the second JET bottle.

To do this, unfasten the levers clockwise to remove the trigger.

Before following the assembly instructions above, give the dip tubes a rinse with tap water and dry with a paper towel.

After use, empty JET bottles and caps can be sent for recycling.

### STORAGE CONDITIONS

JET can be kept for 2 years before being used, and should be stored in a cool dry place, out of direct sunlight, between 10°C - 35°C.

Once opened, JET has an in-use life of 6 months. You can track this by noting the day you opened JET on the bottle label.