

# DUEL



DUEL is a lightweight and portable high-level disinfectant, designed for everyday tasks.

DUEL is proven sporicidal in 30 seconds, making it the fastest high-level disinfectant for surfaces available today.

It is designed specifically for the high-level disinfection of high-touch, non-invasive medical device surfaces in patient bed spaces, intensive care units, isolation wards, operating theatres, oncology wards and neonatal wards.



## DUEL IS SPORICIDAL, MYCOBACTERICIDAL, VIRUCIDAL, FUNGICIDAL, AND BACTERICIDAL.

The contact time of DUEL working solution is realistic. Contact time is the time required for a disinfectant to kill microorganisms. Drying time is the time in which a wetted surface dries. A disinfectant is worthless if the contact time is longer than the drying time. DUEL is effective in 30 seconds against microorganisms of concern such as:

- *Clostridium difficile*
- *Mycobacterium terrae*
- Norovirus
- Influenza A virus (H1N1)
- Hepatitis B virus (HBV)
- Human immunodeficiency virus (HIV)
- *Aspergillus brasiliensis*
- Carbapenem-resistant Enterobacteriaceae (CRE)
- *Klebsiella pneumoniae*
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Vancomycin-resistant Enterococci (VRE) *Enterococcus faecium*

## TEST REPORTS: EUROPEAN NORMS (EN)

BACTERIA	EN 13727
	EN 14561
	EN 16615
YEASTS/FUNGI	EN 13624
	EN 14562
	EN 16615
VIRUS	EN 14476
	EN 16615
MYCOBACTERIA	EN 14348
	EN 14563
SPORES	EN 13704

## POWERED BY CHLORINE DIOXIDE

DUEL utilises a proprietary chlorine dioxide (ClO<sub>2</sub>) chemistry, a well-documented and highly effective biocide. Chlorine dioxide is proven effective in preventing biofilm and removing it from surfaces.


DUEL is far superior to surface cleaners and disinfectants based on bleach, chlorine tablets, quaternary ammonium compounds, phenols and aldehydes.

## DUEL IS PORTABLE

One bottle of DUEL weighs only 0.35kg.

It is lightweight and portable, and can be used wherever required.

DUEL is specifically designed for the high-level disinfection of high-touch, non-invasive medical device surfaces in patient bed spaces.



AVAILABLE TO ORDER NOW!

6 x 250ml DUEL PER BOX  
(310 doses/bottle)

To order quote: CCH040101



**cache**<sup>™</sup>  
SURFACE DISINFECTION

For Tristel patent information please visit:  
[www.our-patents.info/tristel](http://www.our-patents.info/tristel)

**Tristel Solutions Ltd,**  
Lynx Business Park, Fordham Road, Snailwell, Cambs, UK, CB8 7NY  
T: +44 (0)1638 721500 - E: [mail@tristel.com](mailto:mail@tristel.com) - W: [www.tristel.com](http://www.tristel.com)

Copyright © Tristel Solutions  
Mkt-Bro-760-1  
March 2019