

# TRIPLE-ENZYMATIC CLEANING FOAM FOR NON-LUMENED INVASIVE AND NON-INVASIVE MEDICAL DEVICES

including:

- skin surface ultrasound probes
- endocavity ultrasound probes
- non-channel endoscopes

Disinfection is only effective  
when a device has undergone  
thorough cleaning.

- ✓ **Triple Enzymatic**  
Meets cleaning standards
- ✓ **Compatible**  
Approved by major manufacturers
- ✓ **Practical**  
No soaking required
- ✓ **Mobile**  
Usable at point of care
- ✓ **Foam**  
Non-aerosol application

Tristel™  
**CLEAN**



# UNRIVALLED COMPATIBILITY

The Tristel CLEAN formula has been used for over twenty years in the Tristel Pre-Clean Wipe. Tristel cleaning chemistry has been approved by most major manufacturers for use on their devices.

- Advanced Tactile Imaging
- Arthrex
- BD (Bard Access)
- BK Medical
- Butterfly Network
- Canon Medical Systems
- Carestream
- Clarius Mobile Health
- DWL Compumedics
- Echosens
- Esaote
- Exact Imaging
- FUJIFILM Healthcare
- FUJIFILM Sonosite
- GE Healthcare
- Koelis
- Mcube
- Medcom
- Mindray
- Philips
- Siemens Healthineers
- Samsung Medical



## HIGH-LEVEL DISINFECTION WITH TRISTEL DUO FOAMS

As a partner product for Tristel DUO foams and DUO WIPES, Tristel CLEAN is designed for use as part of a comprehensive reprocessing procedure.

After removing soil and organic matter from the device with Tristel CLEAN, Tristel DUO disinfectants are used for high-level disinfection.

Tristel DUO foam is a chlorine dioxide formulation which provides sporicidal, mycobactericidal, virucidal, fungicidal, yeasticidal and bactericidal disinfection in only 30 seconds.

### Product Code:

TSL023301 Pack of 6 x 250ml

TSL024501 Pack of 2 x 250ml

Tristel CLEAN is classified as a Class I Medical Device according to UKCA and EU MDR.

# Tristel™

[www.tristel.com/asia-en](http://www.tristel.com/asia-en)

### Hong Kong

Tristel Asia Ltd

T +852 2895 6968

E [customerservicehk@tristel.com](mailto:customerservicehk@tristel.com)

### India

Tristel India Pvt Ltd

T +91 80417 19452

E [mail-india@tristel.com](mailto:mail-india@tristel.com)

### Singapore

Tristel Pte Ltd

T +65 6518 9918

E [mail-sg@tristel.com](mailto:mail-sg@tristel.com)



Join  
us on  
LinkedIn