

# Tristel™ Sporicidal Wipe

## ACHIEVE HIGH-LEVEL DISINFECTION IN 30 SECONDS WITH THE TRISTEL SPORICIDAL WIPES

For non-lumened medical devices, such as:

- Nasendoscopes
- Transvaginal and transrectal ultrasound probes
- Transoesophageal echocardiogram (TOE/TEE) probes
- Laryngoscopes
- Manometry catheters
- Ophthalmic devices

**DESTROYS  
HPV TYPES  
16 & 18 IN  
30 SECS\***



### **FAST**

Provides high-level disinfection in 30 seconds, increasing patient throughput.



### **MOBILE**

For use at the point-of-care, eliminating the need to transport to central sterilisation and saving valuable time.




### **ECONOMICAL**

No investment required in machines and maintenance contracts.

## TRISTEL CHLORINE DIOXIDE (ClO<sub>2</sub>)

Tristel's proprietary chlorine dioxide (ClO<sub>2</sub>) chemistry is trusted globally in healthcare settings for its fast-acting, easy-to-use and effective disinfection across diverse medical fields.

ClO<sub>2</sub> eliminates pathogens through electron exchange, stealing electrons from the microorganism's structures. Due to this reaction mechanism, microorganisms cannot develop resistance to chlorine dioxide.



### PROVEN EFFECTIVE IN 30 SECONDS AGAINST:

Bacteria	EN 13727
	EN 14561
	EN 16615
Yeasts/Fungi	EN 13624
	EN 14562
	EN 16615
Viruses	EN 14476
	EN 17111
Mycobacteria	EN 14348
	EN 14563
Spores	EN 17126
	EN 17846

## COMPATIBLE

With instruments from major manufacturers, including but not limited to:

- Arthrex
- Canon Medical Systems
- Philips
- ATOMS MedizinTechnik
- GE Healthcare
- Samsung Healthcare
- BD (Bard Access, CareFusion)
- Karl Storz
- Siemens Healthineers
- BK Medical
- Mindray
- Verathon
- Pentax


## EASY-TO-USE



ACTIVATE



WIPE



WAIT

## PRODUCT OPTIONS

### TRISTEL SPORICIDAL WIPES

- 1 box contains 50 Wipes and 1 bottle of Activator Foam

*Note: The Tristel Activator Foam is part of the Tristel Sporcidal Wipes and is not available separately.*



### TRISTEL TRIO WIPES SYSTEM

- 1 box contains 50 decontamination procedures, available with or without the Quality Audit Trail Record Book.

Tristel 3T is available separately.



More information about the Tristel Trio Wipes System, such as safety data sheets, reports, publications on materials and studies, is available on request or online at [www.tristel.com/asia-en](http://www.tristel.com/asia-en)

**Tristel**™

**Hong Kong:** Tristel Asia Limited  
T +852 2895 6968 E [customerservicehk@tristel.com](mailto:customerservicehk@tristel.com)  
**Malaysia:** Tristel Malaysia Sdn. Bhd.  
T +60 37 931 7791 E [mail-my@tristel.com](mailto:mail-my@tristel.com)  
W [www.tristel.com/asia-en](http://www.tristel.com/asia-en)

**India:** Tristel India Pvt Ltd  
T +91 80417 19452 E [mail-india@tristel.com](mailto:mail-india@tristel.com)  
**Singapore:** Tristel Pte Ltd  
T +65 6518 9918 E [mail-sg@tristel.com](mailto:mail-sg@tristel.com)

Join  
us on  
LinkedIn

