



ORGANISM	TEST METHOD	TEST TYPE	CONDITIONS
SPORICIDAL			
<i>Bacillus subtilis</i>	EN 13704	Suspension	Clean 1
<i>Bacillus cereus</i>			
<i>Clostridioides difficile</i>			Clean 1 & Dirty 1

MYCOBACTERICIDAL/ TUBERCULOCIDAL			
<i>Mycobacterium terrae</i>	EN 14348	Suspension	Clean 1
<i>Mycobacterium avium</i>			

VIRUCIDAL			
Poliovirus Type 1	EN 14476	Suspension	Clean 1 & Dirty 1
Adenovirus Type 5			
Norovirus			
Coronavirus (SARS-CoV-2)*			Dirty 1

FUNGICIDAL/ YEASTICIDAL			
<i>Aspergillus brasiliensis</i>	EN 13624	Suspension	Clean 1
<i>Candida albicans</i>			

BACTERICIDAL			
<i>Staphylococcus aureus</i>	EN 14561	Carrier	Clean 1
<i>Enterococcus hirae</i>			
<i>Pseudomonas aeruginosa</i>			
<i>Staphylococcus aureus</i>	EN 13727	Suspension	Clean 1 & Dirty 1
<i>Enterococcus hirae</i>			
<i>Pseudomonas aeruginosa</i>			
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)			

CONDITIONS KEY:

Clean 1: 0.3 g/l Bovine albumin

Dirty 1: 3g/l Bovine albumin + 3ml/l Blood erythrocytes

* A representative sample of Tristel chlorine dioxide chemistry has been tested following EN 14476:2013+A2:2019 at a concentration of 20 parts per million (ppm). Tristel Fuse for Surfaces has a chlorine dioxide concentration greater than 20 ppm at the point of use.