

TRIPLE- ENZYMATIC CLEANING FOAM FOR MEDICAL DEVICES

Tristel™
CLEAN

SAME
PRODUCT,
NEW
LOOK

Suitable for a standalone cleaning process, or prior to disinfection

- ✓ **Triple Enzymatic**
Targeting all major organic soils
- ✓ **Compatible**
Approved by major manufacturers
- ✓ **Practical**
No soaking required
- ✓ **Mobile**
Usable at point of care
- ✓ **Foam**
Non-aerosol application



Endocavity
Probes



Skin Surface
Transducers



Manometry
Catheters



Ophthalmic
Devices



Laryngoscopes



Non-channelled
endoscopes

UNRIVALLED COMPATIBILITY

The Tristel CLEAN triple enzymatic formula has been used for over twenty years in the Tristel Pre-Clean Wipe. Triple enzymatic formulations outperform single and dual enzymes products because they target all major organic soils – **Protease for proteins, Lipase for fats and Amylase for carbohydrates.** Tristel cleaning chemistry has been approved by most major device manufacturers for use on their devices.

APPLICATIONS

Tristel CLEAN can be used as a standalone cleaning product, or prior to disinfection with Tristel DUO Foams, Tristel Trio Wipes System or Tristel Fuse for Stella.



THE PERFECT PARTNERS

Disinfection is only effective when a device has undergone thorough cleaning.

Tristel DUO Foams	Tristel Trio Wipes System	Tristel Fuse for Stella
Designed for use as part of a comprehensive reprocessing procedure. Tristel DUO Foams can be used for disinfection, after removing soil and organic matter from the probe with Tristel CLEAN.	Following patient examination, and prior to transport for decontamination with Trio, Tristel CLEAN can be used for an immediate clean of the medical device.	Tristel CLEAN can be used as a manual and immediate clean for non-lumened medical devices, prior to high level disinfection with Tristel Fuse for Stella chemistry paired with the Stella System.

WHY USE TRISTEL DUO WIPES?

Our dry wipes are designed specifically for the application of CLEAN and DUO foams.



Made of 100% non-woven, low-linting polypropylene



Less absorbent than standard wipes



Does not impede chemical efficacy or scratch surfaces



Australia,
Tristel Pty Ltd, Unit 44, 328 Reserve Road, Cheltenham, VIC 3192, Australia
T 1300 680 898 E mail-au@tristel.com W www.tristel.com/au-en

New Zealand,
23 Birch Avenue, Judea, Tauranga, Bay of Plenty, 3110
T +64 (0) 7 5771560 E mail-nz@tristel.com W www.tristel.com/nz-en

Scan to learn more about Tristel CLEAN!



AU



NZ