

Clean Clever Cuff.

daisy
grip 

daisycuff

Reinventing the traditional blood pressure cuff
with modern infection control practices.

Effortless

fast opening

Reprocess

easily with a
disinfection wipe

Made from skin
friendly material

Attach to
your Digital
Monitoring
System



Designed by doctors

daisycuff solves the infection control problems associated with conventional blood pressure cuffs, while maintaining functionality and ease of use.



why daisycuff?



Materials Matter.

- Manufactured from innovative materials developed and independently inspected as suitable for everyday hospital use.
- Skin-friendly and disinfectant resistant material, enabling fast and complete disinfection between uses, without damaging the cuff.



No Velcro. No Worries.

- The closure mechanism utilises a unique ridged plate and loop system, enabling all functional surfaces to be disinfected.
- Excess material is not required to fold back on itself, saving on unnecessary fabric and surfaces prone to be an infection risk.



Adaptable.

- Sizes range from S (9–16cm) through to XL (32–44cm).
- Adaptors vary to fit most blood pressure machines and manual manometers.



Sustainable.


- Validated for 10,000 uses.
- Durable, reusable and sustainable alternative to single-use devices.

The dirty secrets of blood pressure cuffs

Blood pressure cuffs have long been considered potential vectors for infection. Healthcare-associated infections affect an estimated 170,574 adults in Australian public hospitals each year and cause 7,583 deaths.¹

Published evidence has highlighted contamination rates of blood pressure cuffs, and their potential for infections.

- One study found that 45% of inner sides and 23% of outer sides of cuffs had 100 or more colony-forming units (CFUs) per 25 cm², with the highest rates on the inner side of cuffs from ICUs (83%) and nurses' trolleys (77%). Potentially pathogenic organisms were found on 13% of contaminated cuffs, and genetically linked strains were found on dedicated cuffs used on isolated patients.²
- Another found viable organisms on 100% of cuffs sampled, 58% of them potential pathogens, including MRSA (8%) and *C. difficile* (33%).³
- In Australia, a 2013 study on conventional blood pressure cuffs found high colonisation rates on the inner and outer surfaces of the cuffs in Operating Theatres (76%/86%), the High Dependency Unit (96%/98%) and Emergency (100%/100%). MRSA was found on 2% of all swabs.⁴



References:

1. Lydeamore, M. J., Mitchell, B. G., Bucknall, T., Cheng, A. C., Russo, P. L., & Stewardson, A. J. (2022). Burden of five healthcare associated infections in Australia. *Antimicrobial Resistance & Infection Control*, 11, Article 69. <https://doi.org/10.1086/s13756-022-01109-8>
2. De Gialluly, C., Morange, V., de Gialluly, E., Loulergue, J., van der Mee, N., & Quentin, R. (2006). Blood pressure cuff as a potential vector of pathogenic microorganisms: A prospective study in a teaching hospital. *Infection Control & Hospital Epidemiology*, 27(9), 940–943. <https://doi.org/10.1086/507284>
3. Walker, N., Gupta, R., & Cheesbrough, J. (2006). Blood pressure cuffs: Friend or foe? *Journal of Hospital Infection*, 63(2), 167–169. <https://doi.org/10.1016/j.jhin.2005.10.019>
4. Grewal, H., Varshney, K., Thomas, L. C., Kok, J., & Shetty, A. (2013). Blood pressure cuffs as a vector for transmission of multi-resistant organisms: Colonisation rates and effects of disinfection. *Emergency Medicine Australasia*, 25(3), 222–226. <https://doi.org/10.1111/1742-6723.12076>

Distributed by:

Tristel™

Australian Sponsor:

Tristel Pty Ltd, Unit 44, 328 Reserve Road,
Cheltenham, VIC 3192, Australia
T 1300 680 898 E mail-au@tristel.com
W www.tristel.com/au-en/

Tristel New Zealand:

Tristel New Zealand Ltd, 23 Birch Avenue,
Judea, Tauranga 3110, New Zealand
T +64 (0) 7 5771560 E mail-nz@tristel.com
W www.tristel.com/nz-en/

Scan to
find out
more

