

TRISTEL DUO OPH MATERIAL COMPATIBILITY WITH OPHTHALMIC MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model
Clarity Medical System (Natus)	Ophthalmic Imaging Lenses	RetCam	No Models Specified
DGH Technologies	A-Scan Probe	Scanmate A	DGH 6000
	B-Scan Probe	ScanMate B	DGH 8000
	Desktop Pachymeter Probe	Pachmate	DGH 55, DGH 55B
	Handheld Pachymeter Probe		DGH 55, DGH 55B
Haag-Streit	Contact Glasses	No Series Specified	630 L, 930 L, 906 L, 907 L, Retina 145 L, CGAL, CGIL, CGPL, CGRL, CGVL, 1210, 630, 630 S, 901, 901 S, 902, 902 S, 903, 903 S, 904, 904 S, 905, 905 S, 906, 906 S, 907, 907 S
	Tonometer		All Measuring Prisms
Quantel Medical	Ophthalmic Ultrasound Probes	No Series Specified	No Models Specified
Reichert / Ametek	Ophthalmic Pachymeter	No Series Specified	iPak Pachymeter
Takagi	Tonometer Prism	No Series Specified	No Models Specified
Volk	Ophthalmic Lens	BIO Lenses	Black & All Colors
			ACS
		Classic Series	Black & All Colors
		Super Series	Black & All Colors
		Digital Series Lenses	Black & All Colors
		Mirrored Lenses	3-Mirror Lenses
			Mini 4-Mirror Lens
			SLT
		G-Series Gonio Lenses	No Models Specified
		Contact Lenses	No Models Specified
		Vitrectomy Surgical Lenses	Traditional
ACS			

Disclaimer: Material compatibility testing performed via different methods including wiping and soaking. Tested materials subjected to varying number of cycles.

For more information contact Tristel team on cooperation@tristel.com