

ORGANISME	TESTMETHODE	TESTTYPE	CONDITIES
SPORICIDE			
Bacillus subtilis	EN 17126	Suspensie	Schoon 1
Bacillus cereus			
MYCOBACTERICIDE/TUBERCULOCIDE			
Mycobacterium terrae	EN 14563	Carrier	Schoon 1, Vuil 1 & Vuil 2
Mycobacterium avium			Schoon 1 & Vuil 1
Mycobacterium terrae	EN 14348	Suspensie	Schoon 1 & Vuil 2
Mycobacterium avium			Schoon 1
VIRUCIDE			
Adenovirus Type 5	EN 17111	Carrier	Schoon 1 & Vuil 1
Murine Norovirus			
Poliovirus Type 1	EN 14476	Suspensie	Schoon 1 & Vuil 1
Adenovirus Type 5			
Murine Norovirus			
FUNGICIDE/YEASTICIDE			
Candida albicans	EN 14562	Carrier	Schoon 1 & Vuil 2
Aspergillus brasiliensis	EN 13624	Suspensie	Schoon 1 & Vuil 1
Candida albicans			
BACTERICIDE			
Staphylococcus aureus	EN 14561	Carrier	Schoon 1, Vuil 1 & Vuil 2
Enterococcus hirae			
Pseudomonas aeruginosa			
Staphylococcus aureus	EN 13727	Suspensie	Schoon 1 & Vuil 2
Enterococcus hirae			
Pseudomonas aeruginosa			

Bijkomende Testen

ORGANISME	TESTMETHODE	TESTTYPE	CONDITIES
BACTERIËLE SPOREN			
<i>Clostridium sporogenes</i>	EN 14561/AOAC 966.04	Carrier	Vuil 2
<i>Bacillus subtilis</i>			
MYCOBACTERIËN			
<i>Mycobacterium terrae</i>	DGHM	Carrier	Schoon 1
<i>Mycobacterium avium</i>			
VIRUSSEN			
Bovine Coronavirus	EN 17111	Carrier	Vuil 1
Humaan papillomavirus (HPV) (Polyomavirus SV40 gebruikt als substituu)			
Poliovirus Type 1	ASTM E1053	Suspensie	Vuil 2
Herpes-simplexvirus (HSV) Type 1			
Poliovirus Type 1	DVV/RKI	Suspensie	Schoon 2 & Vuil 3
Adenovirus Type 5			
Murine Norovirus			
Humaan papillomavirus (HPV) (Polyomavirus SV40 gebruikt als substituu)			
Vaccinia Virus			
SCHIMMELS/GISTEN			
<i>Candida albicans</i>	DGHM	Carrier	Schoon 1
<i>Trichophyton interdigitale</i>	AOAC 955.15	Carrier	Vuil 2
<i>Trichophyton interdigitale</i>	AOAC 955.17	Suspensie	Vuil 2

Bijkomende Testen

ORGANISME	TESTMETHODE	TESTTYPE	CONDITIES
CONDITIE			
<i>Klebsiella pneumoniae</i>	EN 14561	Carrier	Schoon 1
Vancomycine-resistente Enterococci (VRE) <i>Enterococcus faecium</i>			
<i>Staphylococcus aureus</i>	DGHM	Carrier	Schoon 1
<i>Enterococcus hirae</i>			
<i>Pseudomonas aeruginosa</i>			
<i>Escherichia coli</i>	EN 13727	Suspensie	Schoon 1

Legende van de condities:

Schoon 1: 0.3g/l Bovine albumine

Schoon 2: Aqua bidest

Vuil 1: 3g/l Bovine albumine + 3ml/l bloed erythrocyten

Vuil 2: 5% Bloedserum

Vuil 3: 10% Bloedserum