

INTRODUCING THE FUTURE OF INSTRUMENT TRANSPORTATION

**TRISTEL PROTECT
BAGS ARE DESIGNED
SPECIFICALLY FOR THE
TRANSPORTATION AND
SHORT-TERM STORAGE
OF INSTRUMENTS, WITH
TRACEABILITY IN MIND**



Tristel™

Protect



72-HOUR STORAGE

WITH TRISTEL PROTECT INSTRUMENTS CAN BE STORED FOR UP TO 72-HOURS* BEFORE USE, REDUCING THE NUMBER OF REPEATED DECONTAMINATION CYCLES.

*after decontamination.



CUTS CAPITAL INVESTMENT

TRISTEL PROTECT REQUIRES MINIMAL CAPITAL INVESTMENT WHEN COMPARED TO A TRADITIONAL TRAY TRANSPORTATION SYSTEM.



SMALL FOOTPRINT

TRISTEL PROTECT'S COMPACT SIZE SAVES SPACE DURING TRANSPORT AND STORAGE.



EASY TO USE

TRISTEL PROTECT BAGS FEATURE AN OPEN-VOID SEAL, 72-HOUR TRACKING LABEL AND TRACEABILITY POCKET, PROVIDING PEACE OF MIND AT EVERY STAGE.



OPEN-VOID SEAL



72-HOUR TRACKING LABEL



TRACEABILITY POCKET

AVAILABLE NOW THROUGH THE NHS SUPPLY CHAIN USING THESE CODES:

SMALL PROTECT BAGS **FAL18777**

LARGE PROTECT BAGS **FAL18356**

MANUFACTURED BY:
TRISTEL SOLUTIONS LIMITED
T 01638 721500
E MAIL@TRISTEL.COM
W WWW.TRISTEL.COM

Tristel[™]
WE HAVE CHEMISTRY.

For Tristel patent information please visit:
<http://www.our-patents.info/tristel>
Copyright © Tristel Solutions
November 2019 - Mkt-Fly-1014-1