

FUSE



2-IN-1 DISINFECTION FOR FLOORS


Designed for use on floors, FUSE combines cleaning action with the sporicidal power of Tristel ClO_2 .

FUSE is proven effective against a wide range of microorganisms, including *Clostridioides difficile* and multi-drug resistant pathogens, in 5 minutes.

One carton of FUSE delivers 200 litres of mopping solution. That is enough to cover 4000m², the surface area of 20 tennis courts!



**BURSTABLE
SACHET
NO MEASURING
NEEDED**



**SPORICIDAL
DISINFECTION
POWERED BY
TRISTEL ClO_2**



**ONE
DILUTION
NO ROOM
FOR ERROR**



**FAST
ACTING
EFFECTIVE IN
5 MINUTES**

POWERED BY TRISTEL ClO₂

FUSE is built upon the same chlorine dioxide (ClO₂) chemistry as the Tristel Medical Device Decontamination products.

ClO₂ is effective even at low concentrations. With alternative mopping products such as chlorine, 25 tablets would need to be diluted into 5 litres of water to reach the same sporicidal efficacy as FUSE. This can result in a solution that is toxic for staff and corrosive for surfaces.

FUSE is delivered easily and rapidly by bursting a sachet and diluting its contents into 5 litres of water. The ClO₂ solution generated is effective at one concentration and in one contact time.

No counting, weighing or calculations needed!




PROVEN EFFECTIVE IN 5 MINUTES, AGAINST:

SPORICIDAL	EN 17126
MYCOBACTERICIDAL	EN 14348
VIRUCIDAL	EN 14476
YEASTICIDAL	EN 13624
	EN 16615
FUNGICIDAL	EN 13624
BACTERICIDAL	EN 13727
	EN 16615

<i>Bacillus cereus</i>	<i>Candida auris</i>
<i>Bacillus subtilis</i>	<i>Trichophyton interdigitale</i>
<i>Clostridioides difficile</i>	Carbapenem-resistant Enterobacteriaceae (CRE)
Adenovirus Type 5	<i>Klebsiella pneumoniae</i>
Murine Norovirus	Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)
Poliovirus Type 1	
SARS-CoV-2*	Multidrug-resistant <i>Acinetobacter baumannii</i> (MDRAB)
Tuberculosis (TB) causing <i>Mycobacteria</i>	Vancomycin-resistant Enterococci (VRE)
<i>Candida albicans</i>	<i>Enterococcus faecium</i>

EASY TO USE BURST, MIX & MOP

Each sachet of FUSE provides you with 5 litres of sporicidal disinfectant solution.
Simply burst, mix and mop!



ALSO AVAILABLE



JET PRO

JET PRO is a sporicidal chlorine dioxide foam which continues to work after application thanks to its unique formulation.



JET

JET is a sporicidal chlorine dioxide foam. Each pull of the trigger targets surfaces where microbes are a risk.



TANK

TANK is a storage and distribution system powered by either Detergent or ClO₂ Sporicidal Disinfectant capsules.



40 Sachet Carton:
CCH080101

ORDER NOW

Our team can offer on-site training, product rollout and audit support for your hospital. Get in touch to find out more.