

**TRISTEL DUO OPH MATERIAL COMPATIBILITY WITH OPHTHALMIC MEDICAL DEVICES**

Manufacturer	Medical Device Group	Medical Device Type	Series	Model
Accutome (Keeler)	Ophthalmic ultrasound probe	Pachymeter	PachPen	n/a
		A-Scan probe	n/a	
		B-Scan probe		
Clarity Medical System (Natus)	Ophthalmic lens / prism	Ophthalmic lens	RetCam	n/a
DGH Technologies	Ophthalmic ultrasound probe	A-Scan probe	n/a	DGH 6000 - Scanmate A
		B-Scan probe		DGH 8000 - Scanmate B
		Pachymeter		DGH 555 - Pachette 3
				DGH 555B - Pachette 4
				DGH 55 - Pachmate
				DGH 55B - Pachmate 2
Ellex	Ophthalmic ultrasound probe	A-Scan probe	Eye Cubed / Eye One	Biometry A-Scan 10 MHz
		B-Scan probe		Diagnostic A-Scan 8 MHz
		A-Scan probe		Anterior B-Scan 40 MHz
				Posterior B-Scan 10 MHz
	B-Scan probe	Eye Prime	Biometry A-Scan 10 MHz, Diagnostic A-Scan 8 MHz	
			Posterior B-Scan 12 MHz, Posterior B-Scan 18 MHz	
	Ophthalmic accessory	Immersion shell	Eye Cubed / Eye One / Eye Prime	n/a
		Scleral shell		n/a
Haag-Streit	Ophthalmic lens / prism	Contact glass	n/a	630 L, 930 L, 906 L, 907 L, Retina 145 L, CGAL, CGIL, CGPL, CGRL, CGVL, 1210, 630, 630 S, 901, 901 S, 902, 902 S, 903, 903 S, 904, 904 S, 905, 905 S, 906, 906 S, 907, 907 S
		Tonometer prism		Measuring prisms for Haag-Streit applanation tonometer
Nidek	Ophthalmic lens / prism	Ophthalmic lens	Gonioscope GS-1	MM prism
	Ophthalmic ultrasound probe	A-Scan probe	n/a	n/a
		B-Scan probe	n/a	n/a
		Pachymeter	n/a	n/a
Quantel Medical	Ophthalmic ultrasound probe	B-Scan probe	n/a	n/a
		Biometry probe		
		Linear probe		
		Pachymeter		
	Ophthalmic accessory	Prager shell		
Reichert / Ametek	Ophthalmic ultrasound probe	Pachymeter	n/a	n/a

**TRISTEL DUO OPH MATERIAL COMPATIBILITY WITH OPHTHALMIC MEDICAL DEVICES**

Manufacturer	Medical Device Group	Medical Device Type	Series	Model
Takagi	Ophthalmic lens / prism	Tonometer prism	n/a	n/a
Tomey	Ophthalmic ultrasound probe	A-Scan probe	n/a	n/a
		B-Scan probe		
		Biometry probe		
		Pachymeter		
	UBM Probe			
	Ophthalmic accessory	Eye Cup		
Volk	Ophthalmic lens / prism	Ophthalmic lens	BIO Lenses	Black & All Colours
				ACS
			Classic Series	Black & All Colours
			Super Series	Black & All Colours
			Digital Series Lenses	Black & All Colours
			Mirrored Lenses	3-Mirror Lenses
				Mini 4-Mirror Lens
				SLT
			G-Series Gonio Lenses	n/a
			Contact Lenses	n/a
	Vitrectomy Surgical Lenses	Traditional		
ACS				
	Ophthalmic accessory	Sterilisation case	Sterilisation Cases	n/a

**Disclaimer:** Material compatibility testing performed via different methods including wiping and soaking. Tested materials subjected to varying number of cycles. Please refer to latest device manufacturer information on compatible cleaners and disinfectants. For more information contact the Tristel team on [devices@tristel.com](mailto:devices@tristel.com)