

Mycobactericidal	EN 14563	Carrier	<i>Mycobacterium avium</i>	MGS	Dirty 1	30 sec.
			<i>Mycobacterium terrae</i>	MGS	Dirty 1	

Virucidal	EN 14476	Suspension	Poliovirus Type 1	Dr Brill & Dr Steinmann	Clean 1	30 sec.
			Adenovirus Type 5	Dr Brill & Dr Steinmann	Clean 1	
			Norovirus (using murine norovirus surrogate)	Dr Brill & Dr Steinmann	Clean 1	
			Influenza (H1N1)	Eurofins	Dirty 4	
	ASTM E-1053	Surface	Poliovirus Type 1	Accuratus Laboratories	Dirty 2	
			Adenovirus Type 5	Accuratus Laboratories	Dirty 2	
			Herpes Simplex Virus Type 1	Accuratus Laboratories	Dirty 2	
	DVV and RKI	Suspension	Poliovirus Type 1	Mikrolab GMBH	Clean 2	
				Mikrolab GMBH	Dirty 3	
			Adenovirus Type 5	Mikrolab GMBH	Clean 2	
				Mikrolab GMBH	Dirty 3	
			Human Papillomavirus (using polyoma virus SV40 surrogate)	Mikrolab GMBH	Clean 2	
Mikrolab GMBH				Dirty 3		
Vaccinia virus			Mikrolab GMBH	Clean 2		
			Mikrolab GMBH	Dirty 3		

Fungicidal/Yeasticidal	EN 14562	Carrier	<i>Candida albicans</i>	Biotech Germande	Clean 1	30 sec.
			<i>Aspergillus brasiliensis</i> (formerly <i>niger</i>)	Hygiene Nord GMBH	Clean 1	
	EN 13624	Suspension	<i>Candida albicans</i>	Hygiene Nord GMBH	Clean 1	
			<i>Aspergillus brasiliensis</i> (formerly <i>niger</i>)	MGS	Clean 1	

Bactericidal	EN 14561	Carrier	<i>Staphylococcus aureus</i>	Biotech Germande	Clean 1	30 sec.
			<i>Pseudomonas aeruginosa</i>	Biotech Germande	Clean 1	
			<i>Enterococcus hirae</i>	Biotech Germande	Clean 1	
			Carbapenemase-resistant Enterobacteriaceae (CRE) <i>Klebsiella pneumoniae</i>	MGS	Clean 1	
			Vancomycin-resistant Enterococci (VRE) <i>Enterococcus faecium</i>	MGS	Clean 1	
	EN 13727	Suspension	<i>Staphylococcus aureus</i>	HIRL	Dirty 4	
			<i>Pseudomonas aeruginosa</i>	HIRL	Dirty 4	
			<i>Enterococcus hirae</i>	HIRL	Dirty 4	

Clean/Dirty Conditions Key

Clean 1: 0.3 g/L bovine albumin
Clean 2: Aqua bidest

Dirty 1: 5% Blood Serum
Dirty 2: 5% Fetal Bovine Serum
Dirty 3: 10% Fetal Calf Serum
Dirty 4: 0.3% albumin + 0.3% sheep erythrocytes



Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

Australian sponsor: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

New Zealand sponsor: Tristel New Zealand Limited, P.O Box 15451, Tauranga, New Zealand
T +64 (0) 7 5771560 - E info@tristel.co.nz

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Copyright © Tristel Solutions
Mkt-Mic-583-2
November 2018