

# STERILE DISINFECTANTS

FOR ENVIRONMENTAL, GLOVED HANDS AND TRANSFER PROCESS DISINFECTION



**Crystel™**

# STERILE DISINFECTANTS

FOR ENVIRONMENTAL, GLOVED HANDS AND TRANSFER PROCESS DISINFECTION

## PLATINUM

### QUATERNARY AMMONIUM COMPOUND/ BIGUANIDE STERILE DISINFECTANT

Contains surface active agents and a blend of quaternary ammonium compounds (QAC) and biguanide for environmental and surface disinfection.

**RECOMMENDED USE:** Platinum is ideal for the disinfection of walls and floors and is suitable for use on most surfaces found within pharmaceutical production facilities.

**DIRECTIONS FOR USE:** All pack sizes must be applied to pre-cleaned surfaces. For professional use only, further use instructions are detailed on the product label.

**CONCENTRATE:** Must be diluted prior to use.

Add one 100ml concentrate to five litres of water to make a 2% solution.

**READY-TO-USE:** Can be used immediately for spraying, wiping and mopping.

**MICROBIOLOGICAL TEST DATA:** EN1276, EN1650, EN13697.

**PRODUCT TESTS CONFIRMED:** Shelf life testing, On-site analytical method, HPLC – High Performance Liquid Chromatography, TOC – Total Organic Compound, Rinsability Graph, Corrosion Studies, COD – Chemical Oxygen Demand, Nitrogen/Phosphorous content.



50 X 100ml - CRY109  
6 X 900ml RTU - CRY110  
2 X 5L RTU - CRY111

## GOLD

### IPA 70% V/V WITH WATER FOR INJECTION

**RECOMMENDED USE:** Gold is ideal for the disinfection of product contact surfaces, spraying in and gloved hands. This ready-to-use solution is suitable for all surfaces found within pharmaceutical production facilities.

**DIRECTIONS FOR USE:** Sterile packs contain ready-to-use 70% v/v IPA with Water for Injection. For professional use only, further use instructions are detailed on the product label.

**MICROBIOLOGICAL TEST DATA:**  
EN1276, EN1650, EN13697.

**PRODUCT TESTS CONFIRMED:**  
Shelf life testing,  
TOC – Total Organic Compound,  
COD – Chemical Oxygen Demand,  
Nitrogen/Phosphorous content.



6 X 1L RTU - CRY101



6 X 1L RTU - CRY104

## SILVER

### DENATURED ETHANOL 70% V/V WITH WATER FOR INJECTION

**RECOMMENDED USE:** Silver is ideal for the disinfection of product contact surfaces, spraying in and gloved hands. This ready-to-use solution is suitable for all surfaces found within pharmaceutical production facilities.

**DIRECTIONS FOR USE:** Sterile packs contain ready-to-use 70% v/v Denatured Ethanol with Water for Injection (Denatured with 5% IPA). For professional use only, further use instructions are detailed on the product label.

**MICROBIOLOGICAL TEST DATA:**  
EN1276, EN1650, EN13697.

**PRODUCT TESTS CONFIRMED:** Shelf life testing, TOC – Total Organic Compound, COD – Chemical Oxygen Demand, Nitrogen/Phosphorous content.

Each Crystel product comes with its own technical file, Certificate of Analysis, Certificates of Sterility and Certificates of Irradiation.

Use biocides safely. Always read the label and product information before use.

#### Brought to you by Tristel in the United Kingdom

Tristel is a manufacturer of products for infection prevention and contamination control, headquartered in Snailwell, Cambs (UK). Please visit [www.tristel.com](http://www.tristel.com) for more information.