

# DISTEL

c

## DOSING SYSTEM



## ELIMINATE SINGLE-USE PLASTIC WIPES

The DISTEL DOSING SYSTEM is intended for the cleaning and disinfection of surfaces of non-invasive medical devices and equipment, and of hard non-porous environmental surfaces.

A 500ml bottle of DISTEL DOSING SYSTEM concentrate generates measured doses of surface disinfectant solution for application via a Doser or Sprayer.

The DISTEL DOSING SYSTEM working solution is effective against coronaviruses in 2 minutes. Efficacy is also proven against bacteria, yeasts, enveloped viruses and nucleic acids (DNA/RNA)\*

\*Please see the microbiological efficacy summary for specific contact times.


Designed to be repurposed and reused, the DISTEL DOSING SYSTEM keeps the environment in mind.

A 500ml bottle of DISTEL DOSING SYSTEM concentrate delivers 10 x 500ml bottles of working solution.

The DISTEL DOSING SYSTEM can be used with any spreader to create freshly made surface disinfectant wipes on the go. Choose a sustainably sourced spreader to reduce single-use plastic waste.


500ML BOTTLE OF CONCENTRATE

1




BOTTLES OF WORKING SOLUTION

10




5ML DOSER

1000 DOSES



1ML SPRAYER

5000 DOSES



USE THE SPREADER OF YOUR CHOICE



MANUFACTURED IN THE UK



REUSABLE BOTTLES



DESTROY CORONAVIRUSES IN 2 MINUTES

# ONE SOLUTION TWO APPLICATORS

The DISTEL DOSING SYSTEM puts you in control.

Choose from two applicators:

- The Doser, for measured doses of liquid
- The Sprayer, for measured doses of liquid spray




**DOSER**

**SPRAYER**

**THE DISTEL DOSING SYSTEM IS INTENDED FOR THE CLEANING AND DISINFECTION OF SURFACES OF NON-INVASIVE MEDICAL DEVICES AND EQUIPMENT, AND OF HARD NON-POROUS ENVIRONMENTAL SURFACES SUCH AS:**




**WORKSTATIONS  
(MONITORS  
& KEYBOARDS)**



**TROLLEYS  
(I.E. SMALL EQUIPMENT  
& RESUSCITATION)**



**HOISTS**



**HANDRAILS**

Training and certification is an essential part of customer service with Cache.  
Ask about our [Online Training Portal](#) for product training at a time and place convenient to you.

**AVAILABLE TO  
ORDER NOW!**



**500ML DISTEL DOSING SYSTEM  
STARTER PACK CCH050401**

1 x 500ml Bottle of Concentrate - 5 x 500ml Dosing  
Bottles - 3 x 5ml Doser  
3 x 1ml Sprayer - 1 x Working Solution Bottle Label



**500ML DISTEL DOSING  
SYSTEM CCH050501**

1 x 500ml Bottle of Concentrate  
1 x Working Solution Bottle Label




**500ML DISTEL DOSING SYSTEM  
DOSING BOTTLES CCH050601**

6 x 500ml Dosing Bottles



**SPRAYER PACK  
CCH030401**

3 x SPRAYER



**DOSER PACK  
CCH030501**

3 x DOSER