

# Tristel

TM

A disinfection cycle is only effective when the medical device has undergone thorough cleaning.

Tristel Clean for Stella efficiently removes organic substances containing proteins, such as blood.

No intermediate rinsing is required if Tristel Fuse for Stella is used for high-level disinfection.



- Tristel Clean for Stella breaks down organic substances and proteins quickly and effectively
- Safe to use
- Two-year shelf life

Tristel Clean for Stella has an integrated 25ml dosing cap, meaning the correct dosage is easy to achieve for each application.

