

Tristel Duo™

OPH

COMPATIBILITY

TRISTEL DUO OPH DEVICE COMPATIBILITY WITH OPHTHALMIC MEDICAL DEVICES

Manufacturer	Medical Device Group	Medical Device Type	Series	Model
DGH Technologies	Ophthalmic ultrasound probe	A-Scan probe	Scanmate	DGH 6000 - Scanmate A
		B-Scan probe		DGH 8000 - Scanmate B
		Pachymeter	Pachette	DGH 555 - Pachette 3
				DGH 555B - Pachette 4
		Handheld pachymeter	Pachmate	DGH 55 - Pachmate
				DGH 55B - Pachmate 2
Ellex	Ophthalmic ultrasound probe	A-Scan probe	Eye Cubed / Eye One	Biometry A-Scan 10 MHz
				Diagnostic A-Scan 8 MHz
		B-Scan probe		Anterior B-Scan 40 MHz
				Posterior B-Scan 10 MHz
		A-Scan probe	Eye Prime	Biometry A-Scan 10 MHz, Diagnostic A-Scan 8 MHz
				Posterior B-Scan 12 MHz, Posterior B-Scan 18 MHz
	B-Scan probe	Eye Cubed / Eye One / Eye Prime	n/a	
				Immersion shell
Ophthalmic accessory	Scleral shell			
	Haag-Streit	Ophthalmic lens / prism	Contact glass	n/a
Tonometer prism			Measuring prisms for Haag-Streit applanation tonometer	
Keeler Accutome		Ophthalmic ultrasound probe	Handheld pachymeter	PachPen (Tip only)
	A-Scan probe		n/a	
	B-Scan probe			
Natus Medical	Ophthalmic lens / prism	Ophthalmic lens	RetCam	n/a
Nidek	Ophthalmic lens / prism	Ophthalmic lens	Gonioscope GS-1	MM prism
	Ophthalmic ultrasound probe	A-Scan probe	n/a	n/a
		B-Scan probe		
		Pachymeter		

TRISTEL DUO OPH DEVICE COMPATIBILITY WITH OPHTHALMIC MEDICAL DEVICES

Manufacturer	Medical Device Group	Medical Device Type	Series	Model
Quantel Medical	Ophthalmic ultrasound probe	B-Scan probe	n/a	n/a
		Biometry probe		
		Linear probe		
		Pachymeter		
	Ophthalmic accessory	Prager Shell		
Reichert Technology	Ophthalmic ultrasound probe	Pachymeter	n/a	n/a
Takagi	Ophthalmic lens / prism	Tonometer prism	n/a	n/a
Tomey	Ophthalmic ultrasound probe	A-Scan probe	n/a	n/a
		B-Scan probe		
		Biometry probe		
	Ophthalmic accessory	Eye cup		
	Ophthalmic ultrasound probe	Pachymeter		
		UBM probe		
Volk	Ophthalmic lens / prism	Ophthalmic lens	BIO Lenses	Black & All Colors
				ACS
			Classic Series	Black & All Colors
			Super Series	
			Digital Series Lenses	
			Mirrored Lenses	3-Mirror Lenses
				Mini 4-Mirror Lens
				SLT
			G-Series Gonio Lenses	n/a
			Contact Lenses	
			Vitreotomy Surgical Lenses	Traditional
				ACS
	Ophthalmic accessory	Sterilisation case	Sterilisation Cases	n/a

Disclaimer: Material compatibility testing performed via different methods including wiping and soaking. Tested materials subjected to varying number of cycles.

For more information contact the Tristel team on devices@tristel.com

Please refer to latest device manufacturer information on compatible cleaners and disinfectants.

Refer to device manufacturers' instructions and care cards for correct handling of devices during reprocessing.



Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Australian sponsor: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

New Zealand: Tristel New Zealand Limited, P.O Box 15451, Tauranga, New Zealand
T +64 (0)7 5771560 - F +64 (0)7 5771567 - E info@tristel.co.nz

Hong Kong & Taiwan: Tristel Asia Limited, 21st Floor, 168 Electric Road, North Point, Hong Kong
T +852 2895 6968 - F +852 2869 4388

Copyright © Tristel Solutions Ltd
Mkt-Val-430-4
November 2022